



ORDER FORM

Pronestander

EILK02

Prone for success with this angle adjustable stander.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		Prov	
PC		Phone	
Email / Fax			

Pronestander

Size 1 Weight Capacity: 132 lbs | Size 2 Weight Capacity: 176 lbs | Size 3 Weight Capacity: 198 lbs

✓ OPTION	DESCRIPTION	RETAIL
Pronestander		\$
Includes Flexible Laterals, PU Knee Supports, Footplate, Tray, Tray Cover and Activity Bowl.		
Must select pelvic band and sandal size below.		
<input type="checkbox"/> 141QE34	Size 1 Pronestander - Meadow	\$2,840
<input type="checkbox"/> 141QE35	Size 2 Pronestander - Mondeo Blue	\$3,018
<input type="checkbox"/> 141QE36	Size 3 Pronestander - Mondeo Blue	\$3,559
Pelvic Band		\$
<input type="checkbox"/> 141ZZ30	Size 1 Pelvic Band - Meadow	\$265
<input type="checkbox"/> 141ZZ31	Size 2 Pelvic Band - Mondeo Blue	\$291
<input type="checkbox"/> 141ZZ32	Size 3 Pelvic Band - Mondeo Blue	\$316
Additional Knee Supports		\$
<input type="checkbox"/> 141XX12	Size 1 Upholstered Knee Cups - Meadow (Pair)	\$182
<input type="checkbox"/> 141XX13	Size 2 Upholstered Knee Cups - Mondeo Blue (Pair)	\$182
<input type="checkbox"/> 141XX14	Size 3 Upholstered Knee Cups - Mondeo Blue (Pair)	\$182
<input type="checkbox"/> 141XX15	Size 4 Upholstered Knee Cups - Mondeo Blue (Pair)	\$182
Foot Positioning		\$
<input type="checkbox"/> 1418A107	Small Sandals with Straps (Pair)	\$273
<input type="checkbox"/> 1418A108	Medium Sandals with Straps (Pair)	\$291
<input type="checkbox"/> 1418A109	Large Sandals with Straps (Pair)	\$315
<input type="checkbox"/> 1418A110	QTY Sandal Raiser 1" (Each) Requires Sandal Selection.	\$107

Technical Specifications - Imperial

Frame Sizing Specification	Size 1	Size 2	Size 3
Age (approximate)	1 - 5 years	4 - 10 years	9 - 18 years
User Height (Min / Max)	31.5" / 45.3"	43.3" / 57"	55" / 73"
Maximum User Weight	132 lbs	176 lbs	198 lbs
Chest Plate Height	18.5" - 27"	27" - 38"	38" - 53"
Base	32" x 24"	42" x 26"	50" x 28"
Activity Tray	16" x 21"	18" x 23"	20" x 25"

Technical Specifications - Metric

Frame Sizing Specification	Size 1	Size 2	Size 3
Age (approximate)	1 - 5 years	4 - 10 years	9 - 18 years
User Height (Min / Max)	800 mm / 1150 mm	1100 mm / 1450 mm	1400 mm / 1850 mm
Maximum User Weight	60 kgs	80 kgs	90 kgs
Chest Plate Height	470 mm - 680 mm	680 mm - 980 mm	970 mm - 1140 mm
Base	820 mm x 600 mm	1070 mm x 650 mm	1270 mm x 700 mm
Activity Tray	400 mm x 525 mm	450 mm x 575 mm	500 mm x 600 mm

Sunrise Medical Canada Inc.
 1000 Creditstone Road, Unit 2
 Concord, ON, L4K 4P8



MK-129853-CAN Rev. 8 © 04/2025

