



ORDER FORM

February 2024

Pronestander

Prone for success with this angle adjustable stander.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		Prov	
PC		Phone	
Email / Fax			

Pronestander

Size 1 Weight Capacity: 132 lbs | Size 2 Weight Capacity: 176 lbs | Size 3 Weight Capacity: 198 lbs

✓ OPTION	DESCRIPTION	RETAIL
Pronestander		\$
Includes Flexible Laterals, PU Knee Supports, Footplate, Tray, Tray Cover and Activity Bowl.		
Must select pelvic band and sandal size below.		
<input type="checkbox"/> LPS/1-MW	Size 1 Pronestander - Meadow	\$2,798
<input type="checkbox"/> LPS/2-MB	Size 2 Pronestander - Mondeo Blue	\$2,973
<input type="checkbox"/> LPS/3-MB	Size 3 Pronestander - Mondeo Blue	\$3,506
Pelvic Band		\$
<input type="checkbox"/> LPS/1/B-MW	Size 1 Pelvic Band - Meadow	\$261
<input type="checkbox"/> LPS/2/B-MB	Size 2 Pelvic Band - Mondeo Blue	\$287
<input type="checkbox"/> LPS/3/B-MB	Size 3 Pelvic Band - Mondeo Blue	\$311
Additional Knee Supports		\$
<input type="checkbox"/> LPS/1/K-MW	Size 1 Upholstered Knee Cups - Meadow (Pair)	\$179
<input type="checkbox"/> LPS/2/K-MB	Size 2 Upholstered Knee Cups - Mondeo Blue (Pair)	\$179
<input type="checkbox"/> LPS/3/K-MB	Size 3 Upholstered Knee Cups - Mondeo Blue (Pair)	\$179
<input type="checkbox"/> LPS/4/K-MB	Size 4 Upholstered Knee Cups - Mondeo Blue (Pair)	\$179
Foot Positioning		\$
<input type="checkbox"/> 152-1600	Small Sandals with Straps (Pair)	\$269
<input type="checkbox"/> 152-2600	Medium Sandals with Straps (Pair)	\$287
<input type="checkbox"/> 152-3600	Large Sandals with Straps (Pair)	\$310
<input type="checkbox"/> 117-787	QTY Sandal Raiser 1" (Each) Requires Sandal Selection.	\$105

Technical Specifications - Imperial

Frame Sizing Specification	Size 1	Size 2	Size 3
Age (approximate)	1 - 5 years	4 - 10 years	9 - 18 years
User Height (Min / Max)	31.5" / 45.3"	43.3" / 57"	55" / 73"
Maximum User Weight	132 lbs	176 lbs	198 lbs
Chest Plate Height	18.5" - 27"	27" - 38"	38" - 53"
Base	32" x 24"	42" x 26"	50" x 28"
Activity Tray	16" x 21"	18" x 23"	20" x 25"

Technical Specifications - Metric

Frame Sizing Specification	Size 1	Size 2	Size 3
Age (approximate)	1 - 5 years	4 - 10 years	9 - 18 years
User Height (Min / Max)	800 mm / 1150 mm	1100 mm / 1450 mm	1400 mm / 1850 mm
Maximum User Weight	60 kgs	80 kgs	90 kgs
Chest Plate Height	470 mm - 680 mm	680 mm - 980 mm	970 mm - 1140 mm
Base	820 mm x 600 mm	1070 mm x 650 mm	1270 mm x 700 mm
Activity Tray	400 mm x 525 mm	450 mm x 575 mm	500 mm x 600 mm



Sunrise Medical Canada Inc.
237 Romina Drive, Unit 3 • Concord, ON L4K 4V3 • Canada



Customer Service: (800) 263-3390 E-mail: cscanada@sunmed.com www.sunrisemedical.ca

Specifications subject to change without notice. © 2024 Sunrise Medical Canada Inc. MK-129853-CAN Rev 6