

FOUR POWERFUL MOTORS

High-torque motors (one powering each wheel) deliver a top speed of 6 mph to conquer tough terrain and maximize battery life for extended range.

UNBEATABLE CURB CLIMBING

Off-road low pressure knobby tires, articulating front suspension, and four independent motors provide unstoppable climbing power.

SMOOTH, JOLT-FREE RIDE

Exceptional suspension and the low-pressure tires keep all four wheels on the ground to smooth out the impact of bumps.

SEATING OPTIONS

From MPS seating to a traditional rehab seat, a whole range of powered functions are available to support your needs.

PERSONALIZED FOR YOU

Choose from six colors and a wide range of options, including Seat Elevate and accessories like phone holders, chargers, and lights.

ANTI-TIP ROLLER

Anti-tip roller enhances traction and control by minimizing excessive tilting on rough terrain.

LEAVE NO TRACE

Low-pressure tires protect trails, preventing erosion and minimizing environmental impact, unlike heavy tracked chairs that can damage sensitive terrain.

DESIGNED IN AUSTRALIA. ASSEMBLED IN THE USA

Since 1994, Magic Mobility has been designing, engineering, and manufacturing highly robust, all-terrain powerchairs so that adventurous riders can experience a world that's not flat.





TECHNICAL DATA	
DIMENSIONS	40.5" Long x 28" Wide
MAX SPEED	6 mph
MAX RANGE	12 miles
BATTERY	Group 24
DRIVE WHEELS	14" Off-Road
DRIVE WHEEL POSITION	4x4
SEAT WIDTH	13" to 23"
SEAT DEPTH	13" to 24"
SEAT-TO-FLOOR HEIGHT	18.5", 19.5", 20.5"
CHAIR WEIGHT	Starting at 320 lbs.
USER WEIGHT CAPACITY	400 lbs. (With power tilt & elevate: 340 lbs.)
POWER SEAT OPTIONS	50° of CG Tilt/Power Recline/12" Seat Elevator/Power Elevating Legrest
COLORS	Midnight Black, Ocean Blue, Fire Engine Red, Clover Green, Tangerine, Camouflage

Refer to order form as specifications may vary based on configuration.

Learn more















