



Customer				Ship To			
Account				Name			
Account #			Date	Address			
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #		City			Province
Buyer				PC		Phone	
Mark For				Email / Fax			

Please submit orders to: cscanada@sunmed.com

Seating Specifications


Typical WC Weight w/ Batteries	Starting at 320 lbs. Weights vary depending on options
Seat Type	MPS and Rehab Style
Power Seating	50° of CG Tilt / Power Recline / 12" Seat Elevator
Seat Size Range	Width: 13-23" Depth: 13-24"
STF Height Range	17", 18", 19"

Drive Base Specifications

Drive Wheel Position	Mid Wheel Drive	Drive Control	R-Net PM 90A controller
User Weight Capacity	400 lbs. 340 lbs. with Elevate or Tilt	Drive Motors	4-pole 6mph
V6 Base Size: All Terrain	Width 28" Length 41.4" Turning Radius 23.5"	Drive Wheels	14" All Terrain Drive Wheels 14" Hybrid Drive Wheels
V6 Base Size: Hybrid	Width 27.3" Length 41.5" Turning Radius 23.5"	Tires / Casters	Pneumatic or Solid inserts
		Charger	Dual Mode
		Battery Size	Group 24 or 34
		Color Accents	4 Frame Colors
		Ground Clearance	3.5"

Client Measurements

a) Top of head to seat	_____	in
b) Top of shoulder to seat	_____	in
c) Underarm to seat	_____	in
d) Elbow to seat	_____	in
e) Top of leg to seat	_____	in
f) Shoulder width	_____	in
g) Chest width	_____	in
i) Hip width	_____	in
j) Actual sitting depth	_____	in
m) Foot length	_____	in
n) Knee to foot	_____	in
o) Seat width	_____	in



Notes

✓ OPTION	DESCRIPTION			RETAIL	
Seating Options				\$	
<input type="checkbox"/> 600SF4	REHAB Seat and Back - Back Canes and adjustable flat seat pan			NC	
<input type="checkbox"/> 600SF5	MPS Seat and Back - Captain style seating			NC	
<input type="checkbox"/> 600SF6	MPS Back and REHAB Seat - Captain style backrest and adjustable flat seat pan			NC	
Base					
Frontier V6 Base		End Item Number: EIPW60		\$	
<input type="checkbox"/> 600MW67	Magic Mobility - Frontier V6 No Power Seating Options			\$ 23,673.00	
<input type="checkbox"/> 600MW68	Magic Mobility - Frontier V6 Single Power Seating Options			\$ 23,673.00	
<input type="checkbox"/> 600MW69	Magic Mobility - Frontier V6 Multi Power Seating Options			\$ 23,673.00	
<hr/>					
<input type="checkbox"/> 600MW70	Magic Mobility - Frontier V6 No Power Seating Options - HD			\$ 23,673.00	
<input type="checkbox"/> 600MW71	Magic Mobility - Frontier V6 Single Power Seating Options - HD			\$ 23,673.00	
<input type="checkbox"/> 600MW72	Magic Mobility - Frontier V6 Multi Power Seating Options - HD			\$ 23,673.00	
<hr/>					
<input type="checkbox"/> 600MW73	Magic Mobility - Frontier V6 No Power Seating Options - MPS (Captain)			\$ 23,673.00	
<input type="checkbox"/> 600MW74	Magic Mobility - Frontier V6 Single Power Seating Options - MPS (Captain)			\$ 23,673.00	
<input type="checkbox"/> 600MW75	Magic Mobility - Frontier V6 Multi Power Seating Options - MPS (Captain)			\$ 23,673.00	
<hr/>					
<input type="checkbox"/> 600MW76	Magic Mobility - Frontier V6 No Power Seating Options - HD MPS (Captain)			\$ 23,673.00	
<input type="checkbox"/> 600MW77	Magic Mobility - Frontier V6 Single Power Seating Options - HD MPS (Captain)			\$ 23,673.00	
<input type="checkbox"/> 600MW78	Magic Mobility - Frontier V6 Multi Power Seating Options - HD MPS (Captain)			\$ 23,673.00	
Drive Wheels				\$	
<input type="checkbox"/> 600DW80	All Terrain Drive Wheels and Tires			NC	
<input type="checkbox"/> 600DW81	Hybrid Drive Wheels and Tires - Pneumatic			NC	
<input type="checkbox"/> 600DW82	Hybrid Drive Wheels and Tires - Solid			\$ 63.00	
Optional Drive Tire Options				\$	
<input type="checkbox"/> 600A90	Tire Sealant for 2 Pneumatic Drive Tires			\$ 51.00	
Caster Wheel Option				\$	
<input type="checkbox"/> 600CA14	Pneumatic Caster	NC	<input type="checkbox"/> 600CA15	Solid Foam Filled Caster	\$ 131.00
Base Accessories				\$	
<input type="checkbox"/> 600A91	Jack for Drive Tire			\$ 272.00	
Frame Color				\$	
<input type="checkbox"/> 600C2	Black Frame Color	NC	<input type="checkbox"/> 600C21	Mossy Green Frame Color	NC
<input type="checkbox"/> 600C3	Red Frame Color	NC	<input type="checkbox"/> 600C195	Royal Blue Frame Color	NC
Seat To Floor Height				\$	
V6 Seat to Floor Heights start at 17" +/- 0.25"					
<input type="checkbox"/> 600HT44	Lowest Seat To Floor Height			NC	
<input type="checkbox"/> 600HT45	Seat To Floor Height +1"			NC	
<input type="checkbox"/> 600HT46	Seat To Floor Height +2"			NC	
Batteries				\$	
<input type="checkbox"/> 600AG25	Group 24 Gel Battery (Pair) - Installed			\$ 1,202.00	
<input type="checkbox"/> 600AG27	Group 34 Gel Battery (Pair) - Installed			\$ 1,335.00	
Battery Charger				\$	
<input type="checkbox"/> 600BJ28	10A Battery Charger			NC	

Seating: Seat Functions

✓ OPTION	DESCRIPTION	RETAIL
No Power Seat Options		
		\$
<input type="checkbox"/> 600ST86	No Power Seating System - 400 lb. weight capacity	\$ 8,698.00
<input type="checkbox"/> 600ST87	No Power Seating System with Power Functionality - 400 lb. weight capacity Required for lights, power recline, power legrest, steering lock, or future power seat functions	\$ 9,328.00

Power Tilt

Utilized without Power Recline.

<input type="checkbox"/> 600PA64	50° Power Tilt in Space - 340 lb. weight capacity	\$ 10,236.00
----------------------------------	---	--------------

Power Seat Elevate

Only available when Power Tilt or Power Tilt and Recline is selected.

<input type="checkbox"/> 600PI16	12" Adjustable Seat Elevator - 340 lb. weight capacity	\$ 8,917.00
----------------------------------	--	-------------

Power Recline

Utilized without Power Tilt. Requires No Power Seating with Power Functionality.

<input type="checkbox"/> 600PA65	150° Power Recline Only - 300 lb. weight capacity	\$ 7,286.00
----------------------------------	---	-------------

Power Tilt and Recline

<input type="checkbox"/> 600PA66	50° Tilt & 150° Power Recline - 300 lb. weight capacity	\$ 16,165.00
----------------------------------	---	--------------

Seating: Seat Sizes

MPS (Captain) Seat Cushion (width x depth)

MPS Backrest can be paired with the REHAB flat pan. If flat pan is desired, omit the MPS cushion and select the desired REHAB flat pan.

<input type="checkbox"/> 600SI24	16" x 18" Vinyl Seat	NC
<input type="checkbox"/> 600SI25	18" x 18" Vinyl Seat	NC
<input type="checkbox"/> 600SI26	18" x 20" Vinyl Seat	NC
<input type="checkbox"/> 600SI27	20" x 18" Vinyl Seat	NC
<input type="checkbox"/> 600SI28	20" x 20" Vinyl Seat	NC
<input type="checkbox"/> 600SI29	22" x 20" Vinyl Seat	\$ 617.00

REHAB Flat Pan Seat Sizes - Adjustable in width and depth for most sizes

REHAB Seat Pan

Standard seat pan available 13-23" wide. Deep seat pan available 13"-23" wide.

<input type="checkbox"/> 600PN13	Standard Adjustable Flat Seat Pan - 13" to 19" deep	NC
<input type="checkbox"/> 600PN14	Deep Adjustable Flat Seat Pan - 20" to 22" deep	\$ 621.00
<input type="checkbox"/> 600PN15	Deep Adjustable Flat Seat Pan - 23" to 24" deep	\$ 621.00

REHAB Seat Width

<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"	<input type="checkbox"/> 20"	<input type="checkbox"/> 21"	<input type="checkbox"/> 22"	<input type="checkbox"/> 23"
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

REHAB Seat Depth

<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"	<input type="checkbox"/> 20"	<input type="checkbox"/> 21"	<input type="checkbox"/> 22"	<input type="checkbox"/> 23"	<input type="checkbox"/> 24"
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

REHAB Seating Cushion Options

Note cushion size will be determined by Rehab Seat size selected above

<input type="checkbox"/> 600CU59	JAY® Ion® (Sizes 16" - 21" Wide)	\$ 657.00
<input type="checkbox"/> 600CU60	JAY® Ion® (Size 22" Wide)	\$ 698.00
<input type="checkbox"/> 600CU61	JAY® Union® (Sizes 16" - 21" Wide)	\$ 764.00
<input type="checkbox"/> 600CU62	JAY® Union® (Size 22" Wide)	\$ 806.00
<input type="checkbox"/> 600CU64	JAY® Fusion® (Sizes 16" - 21" Wide)	\$ 1,045.00
<input type="checkbox"/> 600CU65	JAY® Fusion® (Size 22" Wide)	\$ 1,188.00

Seating: Seat Sizes

✓ OPTION	DESCRIPTION	RETAIL
Support Pads and Hardware (must pick pad and mount)		\$

Lateral Thigh Support Pads

<input type="checkbox"/> 600KF22	Lateral Thigh Support Pad Short - Right (5x5)	\$ 75.00
<input type="checkbox"/> 600KF23	Lateral Thigh Support Pad Short - Left (5x5)	\$ 75.00

Lateral Thigh Support Mounts

<input type="checkbox"/> 600KF26	Lateral Thigh Removable Mount - Right	\$ 325.00
<input type="checkbox"/> 600KF27	Lateral Thigh Removable Mount - Left	\$ 325.00

Seating: Seatback Options

MPS Backrest Options - Comes standard with Manual Recline	\$
---	----

Backrest Width

<input type="checkbox"/> 600BN46	Standard Backrest Width - 16-20" Wide	NC
<input type="checkbox"/> 600BN47	Wide Backrest Width - 22-23" Wide	\$ 452.00

Backrest Height

<input type="checkbox"/> 600B80	Standard Backrest Height - 22" Tall	NC
<input type="checkbox"/> 600B81	Tall Backrest Height - 26" Tall	\$ 744.00

Backrest Material

<input type="checkbox"/> 600BU106	Cloth Backrest	NC
<input type="checkbox"/> 600BU107	Vinyl Backrest	\$ 410.00

MPS Back Accessories	\$
----------------------	----

<input type="checkbox"/> 600PH33	Push Handle - Needed for O2 holder	\$ 211.00
<input type="checkbox"/> 6009A13	Accessory Mount - Needed for aftermarket headrest	\$ 159.00

REHAB Seating Backrest Options	\$
--------------------------------	----

Backrest will have the width of the MPS seat or flat pan. Please choose your back height. Backrest will come with 2 point hardware with a 500 lb. weight capacity.

JAY J3 Back Contour

<input type="checkbox"/> 600JP01	J3@ PA Back for 16"-20" Wide - 2.2" of lateral thoracic support	\$ 912.00
<input type="checkbox"/> 600JP02	J3@ PA Back for 22" Wide - 2.2" of lateral thoracic support	\$ 1,941.00
<input type="checkbox"/> 600JP11	J3@ PL Back for 16"-20" Wide - 3" of lateral thoracic support	\$ 986.00
<input type="checkbox"/> 600JP13	J3@ PD Back for 16"-20" Wide (18" and 21" shell height only) - 6" of lateral thoracic support	\$ 1,068.00

JAY@ J3@ Back Shell Options and Height Ranges

<input type="checkbox"/> 600JR14	16.5" tall shell (Upper Thoracic Short)	NC
<input type="checkbox"/> 600JR15	18" tall shell (Upper Thoracic Med)	NC
<input type="checkbox"/> 600JR16	19.5" tall shell (Upper Thoracic Tall)	NC
<input type="checkbox"/> 600JR17	21" tall shell (Shoulder Height)	NC
<input type="checkbox"/> 600JR27	22.5" tall shell (Shoulder Height)	NC
<input type="checkbox"/> 600JR28	24" tall shell (Shoulder Height)	NC

REHAB Flip Forward Backrest Option	\$
------------------------------------	----

<input type="checkbox"/> 600OD23	REHAB Flip Forward Backrest Option	\$ 214.00
----------------------------------	------------------------------------	-----------

Seating: Seatback Options			
✓ OPTION		DESCRIPTION	RETAIL
JAY® Lateral Thoracic Swing-Away Supports (Requires Pad and Mount) - J3 back only			\$
JAY® Lateral Thoracic Support Pads			
<input type="checkbox"/>	600JW35	Lateral Thoracic (L x H) 4" x 5" Pad - Right	\$ 75.00
<input type="checkbox"/>	600JW38	Lateral Thoracic (L x H) 4" x 5" Pad - Left	\$ 75.00
<input type="checkbox"/>	600JW37	Lateral Thoracic (L x H) 6" x 4" Pad - Right	\$ 75.00
<input type="checkbox"/>	600JW40	Lateral Thoracic (L x H) 6" x 4" Pad - Left	\$ 75.00
JAY® Lateral Thoracic Swing-Away Mount			
<input type="checkbox"/>	600JW41	Lateral Thoracic Swing Away Bracket - Right	\$ 334.00
<input type="checkbox"/>	600JW42	Lateral Thoracic Swing Away Bracket - Left	\$ 334.00
JAY® Anterior Thoracic Support			\$
<input type="checkbox"/>	600ZR10	Bodypoint Monoflex Medium Chest Strap - Pad 3.5" x 18.5" (total length 49.5")	\$ 227.00
<input type="checkbox"/>	600ZR11	Bodypoint Monoflex Large Chest Strap - Pad 3.5" x 21.5" (total length 52.5")	\$ 227.00
<input type="checkbox"/>	600K10	JAY® Anterior Thoracic Support - Medium ATS8302	\$ 322.00
<input type="checkbox"/>	600K11	JAY® Anterior Thoracic Support - Large ATS8303	\$ 322.00
Pelvic Positioning Options - Select One			\$
Pelvic Positioning			
<input type="checkbox"/>	600A848	Positioning Belt, Auto Style Buckle, 2" x 60"	NC
<input type="checkbox"/>	600A849	Positioning Belt, Aircraft Buckle, Padded 2" x 50"	\$ 103.00
<input type="checkbox"/>	600A850	Bodypoint - Evoflex Pelvic Stabilizer Belt	\$ 159.00
<input type="checkbox"/>	600A851	Bodypoint - Evoflex Pelvic Stabilizer Belt with Extension	\$ 272.00
Headrest Options - MPS Style			\$
<input type="checkbox"/>	600NU12	Cloth Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 552.00
<input type="checkbox"/>	600NU13	Vinyl Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 552.00
Headrest Options - REHAB Style			\$
Headrest Pad			
<input type="checkbox"/>	600HE39	Whitmyer® Plush Headrest Pad, 10"	\$ 377.00
<input type="checkbox"/>	600HE46	Whitmyer® Adjust-A-Plush Headrest Pad Medium (Adjustable from 7 to 16" Wide x 4" Tall)	\$ 445.00
Headrest Pad			
<input type="checkbox"/>	600MH33	Whitmyer® Linx Headrest Mount with Quick Release	\$ 247.00
<input type="checkbox"/>	600MH34	Whitmyer® Linx-2 Headrest Mount with Quick Release	\$ 298.00
Armrest Style Options			\$
<input type="checkbox"/>	600EG53	Flip-Up Removable Armrests (Post Style) - Attaches to base. Height adj 3" & angle adj 90°-155°	\$ 1,788.00
<input type="checkbox"/>	600EG54	Cantilever Armrests (Reclining Flip Back) - Attaches to back. Height adj 4" & angle adj 90°	\$ 2,539.00
<input type="checkbox"/>	600EG55	Cantilever 2 Point Quick Release Post - Attaches front of armrest to seat frame	\$ 3,167.00
Armrest Pad Options			\$
<input type="checkbox"/>	600AB81	10" Waterfall Armrest Pad - Right	NC
<input type="checkbox"/>	600AE81	10" Waterfall Armrest Pad - Left	NC
<input type="checkbox"/>	600AB82	14" Waterfall Armrest Pad - Right	NC
<input type="checkbox"/>	600AE82	14" Waterfall Armrest Pad - Left	NC

Legrest Options			
✓ OPTION		DESCRIPTION	RETAIL
Center Mount Fixed Legrest			\$
Center Mount Fixed Legrest			
<input type="checkbox"/> 600H210		Center Mount (14.5" - 17" Leg Length) and Flip-Up Footplate (11" x 10")	\$ 610.00
Center Mount Fixed Legrest - Accessories			
<input type="checkbox"/> 600DI59		Fixed Center Mount Calf Pads	\$ 697.00
<input type="checkbox"/> 600LS16		Foot Strap - 1.5" wide strap	\$ 120.00
Center Mount Power Elevating Legrest			\$
Center Mount Power Elevating Legrest			
<input type="checkbox"/> 600H213		Center Mount Power Elevating - 2 Piece Footplate - Includes calf pads	\$ 5,583.00
Center Mount Power Elevating Legrest Extension			
<input type="checkbox"/> 600E42		Center Mount Power Elevating Extension - Standard (15" - 17" Leg Length)	NC
<input type="checkbox"/> 600E43		Center Mount Power Elevating Extension - Short (10" - 15" Leg Length)	NC
Center Mount Power Elevating Legrest Footplates			
<input type="checkbox"/> 600F117		Center Mount Power Elevating Footplate - Standard (5" x 12" each)	NC
<input type="checkbox"/> 600F118		Center Mount Power Elevating Footplate - Short (4.5" x 7" each)	NC
Swing-Away Legrest			\$
Swing-Away Legrest			
<input type="checkbox"/> 600H211		Swing-Away Legrest (12" - 17" Leg Length) - flip-up footplate (6" x 8.5")	\$ 877.00
Swing-Away Power Elevating Legrest			\$
Swing-Away Power Elevating Legrest			
<input type="checkbox"/> 600H214		Swing-Away Power Elevating Legrest (12" - 17" Leg Length), flip-up height adj footplate (6" x 8.5")	\$ 6,094.00
S/A Power Legrest Positioning Options			
<input type="checkbox"/> 600HL15		Heel Loops - Attaches to the footplate. Length 11.5-13" and Width 1"	\$ 171.00
Accessories			
<input type="checkbox"/> 6008A93		Waterproof Wheelchair Cover	\$ 222.00
<input type="checkbox"/> 6008A94		Waterproof Joystick Cover	\$ 69.00
<input type="checkbox"/> 6008A95		Medical Necessity Bag - Velcro's to armrest (10.5" x 6" x 1.25")	\$ 51.00
<input type="checkbox"/> 600A221		Deluxe Backpack	\$ 144.00
<input type="checkbox"/> 600A615		Cup Holder - Self-leveling. Attaches to armrest opposite joystick	\$ 287.00
<input type="checkbox"/> 600A719		Articulating Oxygen Holder	\$ 1,492.00
<input type="checkbox"/> 600A612		Cane Holder - Includes push handles	\$ 957.00
<input type="checkbox"/> 6009A15		Luggage Rack for V6 - Attaches to base behind the chair (10.5" x 13")	\$ 499.00
<input type="checkbox"/> 6009A17		Lighting Kit for V6 - 2 headlights and 2 red marker lights with turn signals Not available with flip-down knee supports	\$ 1,777.00
<input type="checkbox"/> 6008A97		Charge Mate - 12 Volt/Car outlet and USB/cell phone charger with max draw of 6 amps	\$ 803.00
<input type="checkbox"/> 600A613		Cell Phone Holder - Adjustable to fit most cell phones	\$ 445.00
<input type="checkbox"/> 600A614		Fishing Rod Holder	\$ 445.00
Be Adaptive Products			
<input type="checkbox"/> 6008A98		Be Adaptive Limited Mobility Shooting Mount - Gun or crossbow	\$ 2,242.00
<input type="checkbox"/> 6008A99		Be Adaptive Trigger Mechanism	\$ 1,242.00
<input type="checkbox"/> 6008A100		Be Adaptive Limited Mobility Fishing Mount	\$ 2,242.00

Electronics											
✓ OPTION		DESCRIPTION						RETAIL			
Joystick										\$	
<input type="checkbox"/>	600J121	Standard LED Joystick (w/o Screen Display)						NC			
<input type="checkbox"/>	600J116	Color Joystick (2.3" Color LCD Screen Display, CJSM1)						NC			
<input type="checkbox"/>	600J118	High Resolution Color Joystick,(CJSM2+BT) 3.5" Color LCD Display						\$	757.00		
Joystick Mount										\$	
<input type="checkbox"/>	600JM34	Fixed on Armrest - Depth adjustable 8.25"						NC			
<input type="checkbox"/>	600JM45	Swing-Away Joystick Mount - Includes latch to lock in the forward position						\$	849.00		
<input type="checkbox"/>	600JM46	Midline Joystick Mount - 2 position mount (Rotation 0° - 180°)						\$	2,050.00		
Joystick Mount - Location										\$	
<input type="checkbox"/>	600HM16	Right Hand Armrest						NC			
<input type="checkbox"/>	600HM17	Left Hand Armrest						NC			
Display Language Selection R-net										\$	
<input type="checkbox"/>	600L1	English	NC	<input type="checkbox"/>	600L3	French	NC	<input type="checkbox"/>	600L7	Spanish	NC
Joystick Handle Type										\$	
Comes standard with PG Cone-Shaped Joystick Handle											
<input type="checkbox"/>	600JH13	Ball Handle						\$	110.00		
<input type="checkbox"/>	600JH26	Bodypoint Grooved Mushroom Joystick Handle						\$	240.00		
<input type="checkbox"/>	600JH24	Bodypoint – Rubber Dome Joystick Handle						\$	256.00		
<input type="checkbox"/>	600JH25	Bodypoint – U-Shaped Joystick Handle (4 inch)						\$	256.00		
Expandable Power Module										\$	
Must select both items											
<input type="checkbox"/>	600EX15	Expandable 90A Controller - R-net						\$	1,099.00		
<input type="checkbox"/>	600EX16	Harness for 90A Expandable Controller						\$	768.00		
Power Seat Function										\$	
<input type="checkbox"/>	600EO37	Thru Drive Control 1 Actuator - Must be selected if ordering 1 actuator function						\$	1,443.00		
<input type="checkbox"/>	600EO38	Thru Drive Control 2+ Actuators - Must be selected if ordering 2 or more actuator functions						\$	2,903.00		
Power Seat Control Option										\$	
<input type="checkbox"/>	600QG11	Ctrl+5 Switch Box with Free Assignable Functions						\$	1,017.00		
<input type="checkbox"/>	600QG12	Ctrl+5 Switch Box with Free Assignable Functions (Multi-power only)						NC			
Ctrl+5 Switch Type										\$	
<input type="checkbox"/>	600EY22	Push Button						NC			
<input type="checkbox"/>	600EY23	Toggle Switches						\$	356.00		
Power Seat Control Option - Ctrl+5 Mounting location										\$	
<input type="checkbox"/>	600EY19	Right Armrest	NC	<input type="checkbox"/>	600EY20	Left Armrest	NC				
Attendant Controls										\$	
<input type="checkbox"/>	6000A05	Dual Control Module for Attendant - Includes accessory mount and push handle						\$	2,460.00		
Attendant Control Mounting Location										\$	
<input type="checkbox"/>	600DN5	Right Back Cane	NC	<input type="checkbox"/>	600DN6	Left Back Cane	NC				
R-net Omni2 - 3.5" LCD Display w/ Built-In Infrared & Bluetooth										\$	
<input type="checkbox"/>	600ON14	R-net Omni2 Only Input Control Module						\$	3,611.00		
<input type="checkbox"/>	600ON15	R-net Omni2 Bundled Input Control Module Available w/ Switch-It Input Devices						NC			
R-net Omni2 -Mounting Location										\$	
<input type="checkbox"/>	600DS23	Right Hand Armrest	NC	<input type="checkbox"/>	600DS24	Left Hand Armrest	NC				
Head Control Packages										\$	
<input type="checkbox"/>	600EU25	Switch-It DualPro® Proportional Head Array - Facial Pads						\$	7,057.00		
<input type="checkbox"/>	600EU26	Switch-It DualPro® Proportional Head Array - Spot Pads						\$	7,057.00		