



Customer				Ship To			
Account				Name			
Account #			Date	Address			
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #		City			Prov
Buyer				PC	Phone		
Mark For				Email / Fax			


Please submit orders to: [cscanada@sunmed.com](mailto:cscanada@sunmed.com)

Seating Specifications			
Power Seating	50° of CG Tilt / Power Recline / 12" Seat Elevator	Typical WC Weight w/ Batteries	Starting at 320 lbs. Weights vary depending on options
Seat Size Range	Width: 13-23"   Depth: 13-24"	Seat Type	MPS and Rehab Style
		STF Height Range	17", 18", 19"

Drive Base Specifications			
Drive Wheel Position	Mid Wheel Drive	Drive Control	R-Net PM 90A controller
User Weight Capacity	400 lbs.   340 lbs. with Elevate or Tilt	Drive Motors	4-pole 6mph
V6 Base Size: All Terrain	Width 28"   Length 41.4"   Turning Radius 23.5"	Drive Wheels	14" All Terrain Drive Wheels   14" Hybrid Drive Wheels
V6 Base Size: Hybrid	Width 27.3"   Length 41.5"   Turning Radius 23.5"	Tires / Casters	Pneumatic or Solid inserts
		Charger	Dual Mode
		Battery Size	Group 24 or 34
		Color Accents	4 Frame Colors
		Ground Clearance	3.5"

Client Measurements

- a) Top of head to seat \_\_\_\_\_ in
- b) Top of shoulder to seat \_\_\_\_\_ in
- c) Underarm to seat \_\_\_\_\_ in
- d) Elbow to seat \_\_\_\_\_ in
- e) Top of leg to seat \_\_\_\_\_ in
- f) Shoulder width \_\_\_\_\_ in
- g) Chest width \_\_\_\_\_ in
- i) Hip width \_\_\_\_\_ in
- j) Actual sitting depth \_\_\_\_\_ in
- m) Foot length \_\_\_\_\_ in
- n) Knee to foot \_\_\_\_\_ in
- o) Seat width \_\_\_\_\_ in



Notes

✓ OPTION		DESCRIPTION	RETAIL		
<b>Seating Options</b>			<b>\$</b>		
<input type="checkbox"/> 600SF4		<b>REHAB</b> Seat and Back - Back Canes and adjustable flat seat pan	NC		
<input type="checkbox"/> 600SF5		<b>MPS</b> Seat and Back - Captain style seating	NC		
<input type="checkbox"/> 600SF6		<b>MPS</b> Back and <b>REHAB</b> Seat - Captain style backrest and adjustable flat seat pan	NC		
<b>Base</b>					
Frontier V6 Base		End Item Number: EIPW60	<b>\$</b>		
<input type="checkbox"/> 600MW67		Magic Mobility - Frontier V6 No Power Seating Options	\$ 23,323.00		
<input type="checkbox"/> 600MW68		Magic Mobility - Frontier V6 Single Power Seating Options	\$ 23,323.00		
<input type="checkbox"/> 600MW69		Magic Mobility - Frontier V6 Multi Power Seating Options	\$ 23,323.00		
<input type="checkbox"/> 600MW70		Magic Mobility - Frontier V6 No Power Seating Options - HD	\$ 23,323.00		
<input type="checkbox"/> 600MW71		Magic Mobility - Frontier V6 Single Power Seating Options - HD	\$ 23,323.00		
<input type="checkbox"/> 600MW72		Magic Mobility - Frontier V6 Multi Power Seating Options - HD	\$ 23,323.00		
<input type="checkbox"/> 600MW73		Magic Mobility - Frontier V6 No Power Seating Options - <b>MPS</b> (Captain)	\$ 23,323.00		
<input type="checkbox"/> 600MW74		Magic Mobility - Frontier V6 Single Power Seating Options - <b>MPS</b> (Captain)	\$ 23,323.00		
<input type="checkbox"/> 600MW75		Magic Mobility - Frontier V6 Multi Power Seating Options - <b>MPS</b> (Captain)	\$ 23,323.00		
<input type="checkbox"/> 600MW76		Magic Mobility - Frontier V6 No Power Seating Options - HD <b>MPS</b> (Captain)	\$ 23,323.00		
<input type="checkbox"/> 600MW77		Magic Mobility - Frontier V6 Single Power Seating Options - HD <b>MPS</b> (Captain)	\$ 23,323.00		
<input type="checkbox"/> 600MW78		Magic Mobility - Frontier V6 Multi Power Seating Options - HD <b>MPS</b> (Captain)	\$ 23,323.00		
<b>Drive Wheels</b>			<b>\$</b>		
<input type="checkbox"/> 600DW80		All Terrain Drive Wheels and Tires	NC		
<input type="checkbox"/> 600DW81		Hybrid Drive Wheels and Tires - Pneumatic	NC		
<input type="checkbox"/> 600DW82		Hybrid Drive Wheels and Tires - Solid	\$ 62.00		
<b>Optional Drive Tire Options</b>			<b>\$</b>		
<input type="checkbox"/> 6008A90		Tire Sealant for 2 Pneumatic Drive Tires	\$ 50.00		
<b>Caster Wheel Option</b>			<b>\$</b>		
<input type="checkbox"/> 600CA14	Pneumatic Caster	NC	<input type="checkbox"/> 600CA15	Solid Foam Filled Caster	\$ 129.00
<b>Base Accessories</b>			<b>\$</b>		
<input type="checkbox"/> 6008A91		Jack for Drive Tire	\$ 268.00		
<b>Frame Color</b>			<b>\$</b>		
<input type="checkbox"/> 600C2	Black Frame Color	NC	<input type="checkbox"/> 600C21	Mossy Green Frame Color	NC
<input type="checkbox"/> 600C3	Red Frame Color	NC	<input type="checkbox"/> 600C195	Royal Blue Frame Color	NC
<b>Seat To Floor Height</b>			<b>\$</b>		
V6 Seat to Floor Heights start at 17" +/- 0.25"					
<input type="checkbox"/> 600HT44		Lowest Seat To Floor Height	NC		
<input type="checkbox"/> 600HT45		Seat To Floor Height +1"	NC		
<input type="checkbox"/> 600HT46		Seat To Floor Height +2"	NC		
<b>Batteries</b>			<b>\$</b>		
<input type="checkbox"/> 600AG25		Group 24 Gel Battery (Pair) - Installed	\$ 1,184.00		
<input type="checkbox"/> 600AG27		Group 34 Gel Battery (Pair) - Installed	\$ 1,315.00		
<b>Battery Charger</b>			<b>\$</b>		
<input type="checkbox"/> 600BJ28		10A Battery Charger	NC		

**Seating: Seat Functions**

✓ OPTION	DESCRIPTION	RETAIL
----------	-------------	--------

**No Power Seat Options** \$

<input type="checkbox"/> 600ST86	No Power Seating System - 400 lb. weight capacity	\$ 8,569.00
<input type="checkbox"/> 600ST87	No Power Seating System with Power Functionality - 400 lb. weight capacity Required for lights, power recline, power legrest, steering lock, or future power seat functions	\$ 9,190.00

**Power Tilt** \$

Utilized without Power Recline.

<input type="checkbox"/> 600PA64	50° Power Tilt in Space - 340 lb. weight capacity	\$ 10,085.00
----------------------------------	---	--------------

**Power Seat Elevate** \$

Only available when Power Tilt or Power Tilt and Recline is selected.

<input type="checkbox"/> 600PI16	12" Adjustable Seat Elevator - 340 lb. weight capacity	\$ 8,785.00
----------------------------------	--	-------------

**Power Recline** \$

Utilized without Power Tilt. Requires No Power Seating with Power Functionality.

<input type="checkbox"/> 600PA65	150° Power Recline Only - 300 lb. weight capacity	\$ 7,178.00
----------------------------------	---	-------------

**Power Tilt and Recline** \$

<input type="checkbox"/> 600PA66	50° Tilt & 150° Power Recline - 300 lb. weight capacity	\$ 15,926.00
----------------------------------	---	--------------

**Seating: Seat Sizes**

**MPS (Captain) Seat Cushion (width x depth)** \$

MPS Backrest can be paired with the REHAB flat pan. If flat pan is desired, omit the MPS cushion and select the desired REHAB flat pan.

<input type="checkbox"/> 600SI24	16" x 18" Vinyl Seat	NC
<input type="checkbox"/> 600SI25	18" x 18" Vinyl Seat	NC
<input type="checkbox"/> 600SI26	18" x 20" Vinyl Seat	NC
<input type="checkbox"/> 600SI27	20" x 18" Vinyl Seat	NC
<input type="checkbox"/> 600SI28	20" x 20" Vinyl Seat	NC
<input type="checkbox"/> 600SI29	22" x 20" Vinyl Seat	\$ 608.00

**REHAB Flat Pan Seat Sizes - Adjustable in width and depth for most sizes** \$

**REHAB Seat Pan**

Standard seat pan available 13"-23" wide. Deep seat pan available 13"-23" wide.

<input type="checkbox"/> 600PN13	Standard Adjustable Flat Seat Pan - 13" to 19" deep	NC
<input type="checkbox"/> 600PN14	Deep Adjustable Flat Seat Pan - 20" to 22" deep	\$ 612.00
<input type="checkbox"/> 600PN15	Deep Adjustable Flat Seat Pan - 23" to 24" deep	\$ 612.00

**REHAB Seat Width**

<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"	<input type="checkbox"/> 20"	<input type="checkbox"/> 21"	<input type="checkbox"/> 22"	<input type="checkbox"/> 23"
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

**REHAB Seat Depth**

<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"	<input type="checkbox"/> 20"	<input type="checkbox"/> 21"	<input type="checkbox"/> 22"	<input type="checkbox"/> 23"	<input type="checkbox"/> 24"
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

**REHAB Seating Cushion Options** \$

Note cushion size will be determined by Rehab Seat size selected above

<input type="checkbox"/> 600CU59	JAY® Ion® (Sizes 16" - 21" Wide)	\$ 647.00
<input type="checkbox"/> 600CU60	JAY® Ion® (Size 22" Wide)	\$ 688.00
<input type="checkbox"/> 600CU61	JAY® Union® (Sizes 16" - 21" Wide)	\$ 753.00
<input type="checkbox"/> 600CU62	JAY® Union® (Size 22" Wide)	\$ 794.00
<input type="checkbox"/> 600CU64	JAY® Fusion® (Sizes 16" - 21" Wide)	\$ 1,030.00
<input type="checkbox"/> 600CU65	JAY® Fusion® (Size 22" Wide)	\$ 1,170.00

**Lateral Thigh Support Pads and Hardware (must pick pad and mount)** \$

**Lateral Thigh Support Pads**

<input type="checkbox"/> 600KF22	Lateral Thigh Support Pad Short - Right (5"x5")	\$ 74.00
----------------------------------	---	----------

<input type="checkbox"/>	600KF23	Lateral Thigh Support Pad Short - Left (5"x5")	\$ 74.00
--------------------------	---------	--	----------

**Lateral Thigh Support Mounts**

<input type="checkbox"/>	600KF26	Lateral Thigh Removable Mount - Right	\$ 320.00
<input type="checkbox"/>	600KF27	Lateral Thigh Removable Mount - Left	\$ 320.00

<input checked="" type="checkbox"/>	OPTION	DESCRIPTION	RETAIL
-------------------------------------	--------	-------------	--------

**Seating: Seatback Options**

MPS Backrest Options - Comes standard with Manual Recline			\$
---	--	--	----

**Backrest Width**

<input type="checkbox"/>	600BN46	Standard Backrest Width - 16-20" Wide	NC
<input type="checkbox"/>	600BN47	Wide Backrest Width - 22-23" Wide	\$ 445.00

**Backrest Height**

<input type="checkbox"/>	600B80	Standard Backrest Height - 22" Tall	NC
<input type="checkbox"/>	600B81	Tall Backrest Height - 26" Tall	\$ 733.00

**Backrest Material**

<input type="checkbox"/>	600BU106	Cloth Backrest	NC
<input type="checkbox"/>	600BU107	Vinyl Backrest	\$ 404.00

MPS Back Accessories			\$
----------------------	--	--	----

<input type="checkbox"/>	600PH33	Push Handle - Needed for O2 holder	\$ 208.00
<input type="checkbox"/>	6009A13	Accessory Mount - Needed for aftermarket headrest	\$ 157.00

REHAB Seating Backrest Options			\$
--------------------------------	--	--	----

Backrest will have the width of the MPS seat or flat pan. Please choose your back height. Backrest will come with 2 point hardware with a 500 lb. weight capacity.

**JAY J3 Back Contour**

<input type="checkbox"/>	600JP01	J3® PA Back for 16"-20" Wide - 2.2" of lateral thoracic support	\$ 899.00
<input type="checkbox"/>	600JP02	J3® PA Back for 22" Wide - 2.2" of lateral thoracic support	\$ 1,912.00
<input type="checkbox"/>	600JP11	J3® PL Back for 16"-20" Wide - 3" of lateral thoracic support	\$ 971.00
<input type="checkbox"/>	600JP13	J3® PD Back for 16"-20" Wide (18" and 21" shell height only) - 6" of lateral thoracic support	\$ 1,052.00

**JAY® J3® Back Shell Options and Height Ranges**

<input type="checkbox"/>	600JR14	16.5" tall shell (Upper Thoracic Short)	NC
<input type="checkbox"/>	600JR15	18" tall shell (Upper Thoracic Med)	NC
<input type="checkbox"/>	600JR16	19.5" tall shell (Upper Thoracic Tall)	NC
<input type="checkbox"/>	600JR17	21" tall shell (Shoulder Height)	NC
<input type="checkbox"/>	600JR27	22.5" tall shell (Shoulder Height)	NC
<input type="checkbox"/>	600JR28	24" tall shell (Shoulder Height)	NC

REHAB Flip Forward Backrest Option			\$
------------------------------------	--	--	----

<input type="checkbox"/>	600OD23	REHAB Flip Forward Backrest Option	\$ 211.00
--------------------------	---------	------------------------------------	-----------

**Seating: Seatback Options**

<input checked="" type="checkbox"/>	OPTION	DESCRIPTION	RETAIL
-------------------------------------	--------	-------------	--------

JAY® Lateral Thoracic Swing-Away Supports (Requires Pad and Mount) - J3 back only			\$
---	--	--	----

**JAY® Lateral Thoracic Support Pads**

<input type="checkbox"/>	600JW35	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Right	\$ 74.00
<input type="checkbox"/>	600JW38	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Left	\$ 74.00
<input type="checkbox"/>	600JW37	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Right	\$ 74.00
<input type="checkbox"/>	600JW40	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Left	\$ 74.00

### JAY® Lateral Thoracic Swing-Away Mount

<input type="checkbox"/>	600JW41	Lateral Thoracic Swing Away Bracket - Right	\$ 329.00
<input type="checkbox"/>	600JW42	Lateral Thoracic Swing Away Bracket - Left	\$ 329.00

### JAY® Anterior Thoracic Support

<input type="checkbox"/>	600ZR10	Bodypoint Monoflex Medium Chest Strap - Pad 3.5" x 18.5" (total length 49.5")	\$ 224.00
<input type="checkbox"/>	600ZR11	Bodypoint Monoflex Large Chest Strap - Pad 3.5" x 21.5" (total length 52.5")	\$ 224.00
<input type="checkbox"/>	600K10	JAY® Anterior Thoracic Support - Medium ATS8302	\$ 317.00
<input type="checkbox"/>	600K11	JAY® Anterior Thoracic Support - Large ATS8303	\$ 317.00

### Pelvic Positioning Options - Select One

#### Pelvic Positioning

<input type="checkbox"/>	600A848	Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	NC
<input type="checkbox"/>	600A849	Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	\$ 101.00
<input type="checkbox"/>	600A850	Bodypoint - Evoflex Pelvic Stabilizer Belt	\$ 157.00
<input type="checkbox"/>	600A851	Bodypoint - Evoflex Pelvic Stabilizer Belt with Extension	\$ 268.00

### Headrest Options - MPS Style

<input type="checkbox"/>	600NU12	Cloth Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 544.00
<input type="checkbox"/>	600NU13	Vinyl Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 544.00

### Headrest Options - REHAB Style

#### Headrest Pad

<input type="checkbox"/>	600HE39	Whitmyer® Plush Headrest Pad 254 mm / 10"	\$ 371.00
<input type="checkbox"/>	600HE46	Whitmyer® Adjust-A-Plush Headrest Pad Medium (Adjustable from 7 to 16" Wide x 4" Tall)	\$ 438.00

#### Headrest Pad

<input type="checkbox"/>	600MH33	Whitmyer® Linx Headrest Mount with Quick Release	\$ 243.00
<input type="checkbox"/>	600MH34	Whitmyer® Linx-2 Headrest Mount with Quick Release	\$ 294.00

### Armrest Style Options

<input type="checkbox"/>	600EG53	Flip-Up Removable Armrests (Post Style) - Attaches to base. Height adj 3" & angle adj 90°-155°	\$ 1,762.00
<input type="checkbox"/>	600EG54	Cantilever Armrests (Reclining Flip Back) - Attaches to back. Height adj 4" & angle adj 90°	\$ 2,501.00
<input type="checkbox"/>	600EG55	Cantilever 2 Point Quick Release Post - Attaches front of armrest to seat frame	\$ 3,120.00

### Armrest Pad Options

<input type="checkbox"/>	600AB80	4" x 16" Tapered Pads - Right	\$ 146.00
<input type="checkbox"/>	600AE80	4" x 16" Tapered Pads - Left	\$ 146.00
<input type="checkbox"/>	600AB81	10" Waterfall Armrest Pad - Right	NC
<input type="checkbox"/>	600AE81	10" Waterfall Armrest Pad - Left	NC
<input type="checkbox"/>	600AB82	14" Waterfall Armrest Pad - Right	NC
<input type="checkbox"/>	600AE82	14" Waterfall Armrest Pad - Left	NC

### Legrest Options

✓ OPTION	DESCRIPTION	RETAIL
----------	-------------	--------

### Center Mount Fixed Legrest

#### Center Mount Fixed Legrest

<input type="checkbox"/>	600H210	Center Mount Flip-Up Footplate	\$ 601.00
--------------------------	---------	--------------------------------	-----------

#### Center Mount Fixed Legrest - Calfpads

<input type="checkbox"/>	600DI59	Fixed Center Mount Calf Pads	\$ 687.00
--------------------------	---------	------------------------------	-----------

**Center Mount Power Elevating Legrest** \$

**Center Mount Power Elevating Legrest**

<input type="checkbox"/> 600H213	Center Mount Power Elevating - 2 Piece Footplate - Includes calf pads	\$ 5,500.00
----------------------------------	---	-------------

**Center Mount Power Elevating Legrest Extension**

<input type="checkbox"/> 600E42	Center Mount Power Elevating Extension - Standard	NC
<input type="checkbox"/> 600E43	Center Mount Power Elevating Extension - Short	NC

**Center Mount Power Elevating Legrest Footplates**

<input type="checkbox"/> 600F117	Center Mount Power Elevating Footplate - Standard	NC
<input type="checkbox"/> 600F118	Center Mount Power Elevating Footplate - Short	NC

**Swing-Away Legrest** \$

**Swing-Away Legrest**

<input type="checkbox"/> 600H211	Side Mount Swing-Away Legrest - 2 piece flip-up footplate (8.5" x 6"). 13" height adjustable	\$ 864.00
----------------------------------	--	-----------

**S/A Non-Power Positioning Options**

<input type="checkbox"/> 600LS15	Heel and Side Support - Attaches to the footplate stem and footplate	\$ 870.00
<input type="checkbox"/> 600LS16	Foot Strap - 1.5" wide strap	\$ 118.00

**Swing-Away Power Elevating Legrest** \$

**Swing-Away Power Elevating Legrest**

<input type="checkbox"/> 600H214	Side Mount Power Elevating Legrest - Swing away legrest with flip-up 2" height adj footplate	\$ 6,004.00
----------------------------------	--	-------------

**S/A Power Legrest Positioning Options**

<input type="checkbox"/> 600HL15	Heel Loops - Attaches to the footplate. Length 11.5-13" and Width 1"	\$ 168.00
----------------------------------	--	-----------

**Accessories** \$

<input type="checkbox"/> 6008A93	Waterproof Wheelchair Cover	\$ 219.00
<input type="checkbox"/> 6008A94	Waterproof Joystick Cover	\$ 68.00
<input type="checkbox"/> 6008A95	Medical Necessity Bag - Velcro's to armrest (10.5" x 6" x 1.25")	\$ 50.00
<input type="checkbox"/> 600A221	Deluxe Backpack	\$ 142.00
<input type="checkbox"/> 600A615	Cup Holder - Self-leveling. Attaches to armrest opposite joystick	\$ 283.00
<input type="checkbox"/> 600A719	Articulating Oxygen Holder	\$ 1,470.00
<input type="checkbox"/> 600A612	Cane Holder - Includes push handles	\$ 943.00
<input type="checkbox"/> 6009A15	Luggage Rack for V6 - Attaches to base behind the chair (10.5" x 13")	\$ 492.00
<input type="checkbox"/> 6009A17	Lighting Kit for V6 - 2 headlights and 2 red marker lights with turn signals Not available with flip-down knee supports	\$ 1,751.00
<input type="checkbox"/> 6008A97	Charge Mate - 12 Volt/Car outlet and USB/cell phone charger with max draw of 6 amps	\$ 791.00
<input type="checkbox"/> 600A613	Cell Phone Holder - Adjustable to fit most cell phones	\$ 438.00
<input type="checkbox"/> 600A614	Fishing Rod Holder	\$ 438.00

**Be Adaptive Products**

<input type="checkbox"/> 6008A98	Be Adaptive Limited Mobility Shooting Mount - Gun or crossbow	\$ 2,209.00
<input type="checkbox"/> 6008A99	Be Adaptive Trigger Mechanism	\$ 1,224.00
<input type="checkbox"/> 6008A100	Be Adaptive Limited Mobility Fishing Mount	\$ 2,209.00

**Electronics**

✓ OPTION	DESCRIPTION	RETAIL
<b>Joystick</b> <span style="float: right;">\$</span>		
<input type="checkbox"/> 600J121	Standard LED Joystick (w/o Screen Display)	NC
<input type="checkbox"/> 600J116	Color Joystick (2.3" Color LCD Screen Display, CJSM1)	NC
<input type="checkbox"/> 600J118	High Resolution Color Joystick,(CJSM2+BT) 3.5" Color LCD Display	\$ 746.00

<b>Joystick Mount</b>					\$						
<input type="checkbox"/>	600JM34		Fixed on Armrest - Depth adjustable 8.25"		NC						
<input type="checkbox"/>	600JM45		Swing-Away Joystick Mount - Includes latch to lock in the forward position		\$ 836.00						
<input type="checkbox"/>	600JM46		Midline Joystick Mount - 2 position mount (Rotation 0° - 180°)		\$ 2,020.00						
<b>Joystick Mount - Location</b>					\$						
<input type="checkbox"/>	600HM16		Right Hand Armrest		NC						
<input type="checkbox"/>	600HM17		Left Hand Armrest		NC						
<b>Display Language Selection R-net</b>					\$						
<input type="checkbox"/>	600L1	English	NC	<input type="checkbox"/>	600L3	French	NC	<input type="checkbox"/>	600L7	Spanish	NC
<b>Joystick Handle Type</b>					\$						
Comes standard with PG Cone-Shaped Joystick Handle											
<input type="checkbox"/>	600JH13		Ball Handle		\$ 108.00						
<input type="checkbox"/>	600JH26		Bodypoint Grooved Mushroom Joystick Handle		\$ 236.00						
<input type="checkbox"/>	600JH24		Bodypoint – Rubber Dome Joystick Handle		\$ 252.00						
<input type="checkbox"/>	600JH25		Bodypoint – U-Shaped Joystick Handle (4 inch)		\$ 252.00						
<b>Expandable Power Module</b>					\$						
Must select both items											
<input type="checkbox"/>	600EX15		Expandable 90A Controller - R-net		\$ 1,083.00						
<input type="checkbox"/>	600EX16		Harness for 90A Expandable Controller		\$ 757.00						
<b>Power Seat Function</b>					\$						
<input type="checkbox"/>	600EO37		THRU DRIVE CNTRL 1 ACTUATOR- Must be selected if ordering 1 actuator function		\$ 1,422.00						
<input type="checkbox"/>	600EO38		THRU DRIVE CNTRL 2+ ACTUATORS- Must be selected if ordering 2 or more actuator functions		\$ 2,860.00						
<b>Power Seat Control Option</b>					\$						
<input type="checkbox"/>	600QG11		Ctrl+5 Switch Box with Free Assignable Functions		\$ 1,002.00						
<input type="checkbox"/>	600QG12		Ctrl+5 Switch Box with Free Assignable Functions (Multi-power only)		NC						
<b>Ctrl+5 Switch Type</b>					\$						
<input type="checkbox"/>	600EY22		Push Button		NC						
<input type="checkbox"/>	600EY23		Toggle Switches		\$ 351.00						
<b>Power Seat Control Option - Mounting location</b>					\$						
<input type="checkbox"/>	600EY19	CTRL+5 MNT ARMREST RIGHT	NC	<input type="checkbox"/>	600EY20	CTRL+5 MNT ARMREST LEFT	NC				
<b>Attendant Controls</b>					\$						
<input type="checkbox"/>	6000A05		Dual Control Module for Attendant - Includes accessory mount and push handle		\$ 2,424.00						
<b>Attendant Control Mounting Location</b>					\$						
<input type="checkbox"/>	600DN5	Right Back Cane	NC	<input type="checkbox"/>	600DN6	Left Back Cane	NC				
<b>R-net Omni-2 - 3.5" LCD Display w/ Built-In Infrared &amp; Bluetooth</b>					\$						
<input type="checkbox"/>	600ON14		R-net Omni-2 Only Input Control Module		\$ 3,558.00						
<input type="checkbox"/>	600ON15		R-net Omni2 Bundled Input Control Module Available w/ Switch-It Input Devices		NC						
<b>R-net Omni-2 -Mounting Location</b>					\$						
<input type="checkbox"/>	600DS23	Right Hand Armrest	NC	<input type="checkbox"/>	600DS24	Left Hand Armrest	NC				
<b>Head Control Packages</b>					\$						
<input type="checkbox"/>	600EU25		Switch-It DualPro® Proportional Head Array - Facial Pads		\$ 6,953.00						
<input type="checkbox"/>	600EU26		Switch-It DualPro® Proportional Head Array - Spot Pads		\$ 6,953.00						