

RGK

TIGA

» Ideal for all aspects of a busy daily life



TIGA - Your customised rigid chair with maximum performance

The TIGA uses advanced technology and high quality material to create the ultimate all-rounder chair.

The TIGA is made from Aluminium 7020 and is perfect to suit your busy lifestyle, made to your individual measurements and personal preferences.

This great performing chair is clean with an attractive and modern design.


Key Features

- » Customized rigid frame Aluminium 7020
- » Available in open and closed frame design
- » Modern look
- » Integrated fix COG or adjustable COG available
- » Integrated fixed backrest or folding back available
- » Ergonomic seat and back optional
- » Wide accessories range
- » Choice of anodised accent colours
- » Airtech tension adjustable seat and back upholstery
- » Max user weight 125 kg

Open Frame




Closed Frame



Choice of anodised accent colours



Integrated lateral supports adjustable. Additional stability for those users that require more lateral support than provided by the upholstery. Helps preventing scoliosis or its progression



Ergonomic seat and back (optional). Increases your stability, functionality and skin protection.



FrontWheel

Enables you to travel easily across a wide range of terrains by keeping your front castor wheels off the floor at all times. Attaches to TIGA, TIGA fx and Hi Lite.