



**Quickie QXi**  
February 2024

EIF6

<b>Mark For:</b> Date: _____ Dealer Acct #: _____ Dealer: _____ Dealer Contact: _____ Dealer Address: _____ Dealer City: _____ Prov: _____ PC: _____ Dealer Phone: ( ) _____ Fax: ( ) _____ Confirmation Email: _____ Confirm Via: <input type="checkbox"/> Fax <input type="checkbox"/> Email	<b>Submitting for:</b> <input type="checkbox"/> Quote <input type="checkbox"/> Order PO#: _____  <p style="text-align: center;"><b>ADDITIONAL SHIPPING INFORMATION</b></p> Ship To: _____ Attention: _____ Address: _____ Address: _____ Ship To City: _____ Prov _____ PC: _____ Ship To Phone: ( ) _____ Fax: ( ) _____
---	---

✓ Option #	Definition	Price	✓ Option #	Definition	Price
------------	------------	-------	------------	------------	-------

**1. FRAME**

Base Model			
<i>Weight capacity 300 lbs</i>			
<input type="checkbox"/>	196WM78	Quickie QXi (Adjustable Axle) For 14" sizes use E1236	\$2,677.00

Frame Style			
<i>See Seat Height Matrix on page 2</i>			
<input type="checkbox"/>	196FS9	Standard Frame	NC
<input type="checkbox"/>	196FS30	Hemi Frame	NC

Transit			
<input type="checkbox"/>	196TO38	Transit System WC19	NC

**2. SIZE / SEATING**

Frame Width			
<input type="checkbox"/>	196W5	14"	NC
<input type="checkbox"/>	196W6	15"	NC
<input type="checkbox"/>	196W7	16"	NC
<input type="checkbox"/>	196W8	17"	NC
<input type="checkbox"/>	196W9	18"	NC
<input type="checkbox"/>	196W10	19"	\$303.00
<input type="checkbox"/>	196W11	20"	\$303.00

Seat Depth			
<i>Frame Depths equals Seat Depths.</i>			
<input type="checkbox"/>	196SD7	14"	NC
<input type="checkbox"/>	196SD8	15"	NC
<input type="checkbox"/>	196SD9	16"	NC
<input type="checkbox"/>	196SD10	17"	NC
<input type="checkbox"/>	196SD11	18"	NC
<input type="checkbox"/>	196SD71	19"	\$303.00
<input type="checkbox"/>	196SD13	20"	\$303.00

Frame Width - Built 4 Me Options			
<i>Not available with EXO standard or Tension adjustable upholstery, Not available with Ballistic Nylon or 3DX Tension adjustable upholstery.</i>			
<input type="checkbox"/>	196W59	B4ME - 21" Width	\$303.00
<input type="checkbox"/>	196W60	B4ME - 22" Width	\$303.00

JAY Cushions (Optional)			
<i>To order any JAY cushions on-chair, please see the JAY price list.</i>			
<i>Cushion name or part number must be filled out with desired width and depth</i>			
<input type="checkbox"/>	196JC20	JAY Bundled Cushion Cushion Name or Part # _____ Width: _____ Depth: _____	Price List

Seat Sling			
<input type="checkbox"/>	196SL1	Standard Seat Sling	NC
<input type="checkbox"/>	196SL6	Omit Seat Sling	NC

Personalized Seating			
<input type="checkbox"/>	196QC01	Surefit	
<input type="checkbox"/>	196QC02	Configurfit	
<input type="checkbox"/>	196QC03	JayFit	



Quickie QXi, Ultra Lightweight Adjustable Folding Wheelchair



Quickie QXi featuring the Jay ION™ cushion

✓	Option #	Definition	Price	✓	Option #	Definition	Price
---	----------	------------	-------	---	----------	------------	-------

**3. CASTERS / SEAT HEIGHT**

Front Seat Height			
<i>See Front Seat Height Matrix. 12.5 &amp; 13.0 only available w/ 60° hanger 12.5" &amp; 13" only available w/ Alum Angle Adj Footplate</i>			
<input type="checkbox"/>	196SH34	12.5"	NC
<input type="checkbox"/>	196SH35	13.0"	NC
<input type="checkbox"/>	196SH40	13.5"	NC
<input type="checkbox"/>	196SH41	14.0"	NC
<input type="checkbox"/>	196SH42	14.5"	NC
<input type="checkbox"/>	196SH12	15.0"	NC
<input type="checkbox"/>	196SH13	15.5"	NC
<input type="checkbox"/>	196SH14	16.0"	NC
<input type="checkbox"/>	196SH15	16.5"	NC
<input type="checkbox"/>	196SH16	17.0"	NC
<input type="checkbox"/>	196SH17	17.5"	NC
<input type="checkbox"/>	196SH18	18.0"	NC
<input type="checkbox"/>	196SH20	18.5"	NC
<input type="checkbox"/>	196SH21	19.0"	NC
<input type="checkbox"/>	196SH23	19.5"	NC

Caster Wheels			
<input type="checkbox"/>	196CW3	4" Low Profile Polyurethane	NC
<input type="checkbox"/>	196CW44	4" x 1.25" Semi Pneumatic	NC
<input type="checkbox"/>	196CW5	5" Low Profile Polyurethane	NC
<input type="checkbox"/>	196CW7	5" x 1.5" Semi Pneumatic	NC
<input type="checkbox"/>	196CW34	6" Polyurethane	NC
<input type="checkbox"/>	196CW35	6" x 1.5" Semi Pneumatic	NC
<input type="checkbox"/>	196CW32	8" Polyurethane	NC
<input type="checkbox"/>	196CW39	8" x 1.5" Semi Pneumatic	\$144.00

Caster Fork			
<input type="checkbox"/>	196CF71	3.5" (standard stem bolt only)	NC
<input type="checkbox"/>	196CF12	5.25"	NC
<input type="checkbox"/>	196CF17	6"	NC
<input type="checkbox"/>	196CF23	7"	NC

Caster Stem Bolt			
<input type="checkbox"/>	196CS1	Standard	NC
<input type="checkbox"/>	196CS5	3/4"	NC
<input type="checkbox"/>	196CS9	1 1/2"	NC

Front Seat Height Matrix								
Front Seat Height Matrix - Standard Stem Bolt								
Caster	3.5" FORK				5 1/4" FORK			
	STD		HEMI		STD		HEMI	
	Min	Max	Min	Max	Min	Max	Min	Max
4 Low Poly	14.5	15.0	12.5	13.0	16.0	17.0	14.0	15.0
4 x 1.25 Semi Pneu	15.0	15.5	13.0	13.5	16.5	17.5	14.5	15.5
5 Low Poly	15.0	15.5	13.0	13.5	16.5	17.5	14.5	15.5
5 x 1.5 Semi Pneu	15.0	15.5	13.0	13.5	16.5	17.5	14.5	15.5
6 Poly	N/A	N/A	N/A	N/A	17.0	18.0	15.0	16.0
6 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	17.0	18.0	15.0	16.0
8 Poly	N/A	N/A	N/A	N/A	19.0	19.0	17.0	17.0
8 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	18.5	19.0	16.5	17.0

Caster	6" FORK				7" FORK			
	STD		HEMI		STD		HEMI	
	Min	Max	Min	Max	Min	Max	Min	Max
4 Low Poly	17.0	18.0	15.0	16.0	18.0	19.0	16.0	17.0
4 x 1.25 Semi Pneu	17.0	18.0	15.0	16.0	18.0	19.0	16.0	17.0
5 Low Poly	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5
5 x 1.5 Semi Pneu	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5
6 Poly	17.5	19.0	15.5	17.0	18.5	19.5	16.5	18.0
6 x 1.5 Semi Pneu	17.5	19.0	15.5	17.0	18.5	19.5	16.5	18.0
8 Poly	19.5	19.5	17.5	18.0	19.5	19.5	17.5	19.0
8 x 1.5 Semi Pneu	18.5	19.5	16.5	17.5	19.0	19.5	17.0	18.5

*Seat squeeze and seat depth can affect FSH availability  
For 3/4" and 1 1/2" Stem bolt matrix see page 5*

**4. REAR WHEELS / REAR SEAT HEIGHTS**

Rear Seat Height			
<i>See Rear Seat Height Matrix</i>			
<input type="checkbox"/>	196RH9	12"	NC
<input type="checkbox"/>	196RH10	12.5"	NC
<input type="checkbox"/>	196RH11	13"	NC
<input type="checkbox"/>	196RH12	13.5"	NC
<input type="checkbox"/>	196RH13	14"	NC
<input type="checkbox"/>	196RH14	14.5"	NC
<input type="checkbox"/>	196RH15	15"	NC
<input type="checkbox"/>	196RH16	15.5"	NC
<input type="checkbox"/>	196RH17	16"	NC
<input type="checkbox"/>	196RH18	16.5"	NC
<input type="checkbox"/>	196RH19	17"	NC
<input type="checkbox"/>	196RH26	17.5"	NC
<input type="checkbox"/>	196RH20	18"	NC
<input type="checkbox"/>	196RH27	18.5"	NC
<input type="checkbox"/>	196RH21	19"	NC

Rear Wheel Size			
<i>See Rear Tire &amp; Wheel Matrix for availability</i>			
<input type="checkbox"/>	196RW3	20" Mag Wheel	NC
<input type="checkbox"/>	196RW8	22" Mag Wheel	NC
<input type="checkbox"/>	196RW1	24" Mag Wheel	NC
<input type="checkbox"/>	196RW103	20" Mag Wheel 5 Spoke	NC
<input type="checkbox"/>	196RW104	22" Mag Wheel 5 Spoke	NC
<input type="checkbox"/>	196RW102	24" Mag Wheel 5 Spoke	NC
<input type="checkbox"/>	196RW90	20" Lite Spoke Wheel	NC
<input type="checkbox"/>	196RW91	22" Lite Spoke Wheel	NC
<input type="checkbox"/>	196RW92	24" Lite Spoke Wheel	NC
<input type="checkbox"/>	196RW93	25" Lite Spoke Wheel	NC
<input type="checkbox"/>	196RW94	26" Lite Spoke Wheel	\$124.00
<input type="checkbox"/>	196RW25	Omit	NC

Rear Seat Height Matrix		
Rear Wheel Size	Std Frame	Hemi Frame
20" Rear Wheel	13.0" - 16.5"	12.0" - 15.0"
22" Rear Wheel	14.0" - 17.5"	12.5" - 16.0"
24" Rear Wheel	15.0" - 18.5"	13.5" - 16.5"
25" Rear Wheel	15.5" - 19.0"	14.0" - 17.0"
26" Rear Wheel	16.0" - 19.5"	14.5" - 18.0"

Rear Wheel Tire			
<i>See Rear Tire &amp; Wheel Matrix for availability</i>			
<input type="checkbox"/>	196RT1	Pneumatic	NC
<input type="checkbox"/>	196RT2	Pneumatic w/ Airless Insert	\$243.00
<input type="checkbox"/>	196RT3	Full Poly	NC
<input type="checkbox"/>	196RT39	High-Pressure Clincher Black	\$103.00
<input type="checkbox"/>	196RT21	Schwalbe Marathon Plus Evolution	\$155.00

Handrim			
<i>Natural Fit only available w/ 22, 24, 25 Lite Spoke. Plastic N/A w/ 25, 26"</i>			
<i>See Appendix for Rear wheel / Handrim matrix</i>			
<input type="checkbox"/>	196HR1	Aluminum Anodized	NC
<input type="checkbox"/>	196HR4	Plastic	\$144.00
<input type="checkbox"/>	196HR22	Natural Fit - Standard Grip	\$368.00
<input type="checkbox"/>	196HR23	Natural Fit - Super Grip	\$428.00
<input type="checkbox"/>	196HR32	Natural Fit - No Thumb Grip	\$549.00
<input type="checkbox"/>	196HR26	The Surge - Oval w/Gription Strip	\$549.00
<input type="checkbox"/>	196HR27	The Surge LT - Oval w/Gription Strip	\$583.00
<input type="checkbox"/>	196HR28	Q-Grip Neoprene Coated	\$464.00

Seat Squeeze Matrix	
Frame Depth	Max Seat Squeeze
14"	2.0"
15" - 16"	2.5"
17" - 18"	2.5"
19" - 20"	3.0"

**4. REAR WHEELS / REAR SEAT HEIGHTS (CONT)**

**Rear Tire & Wheel Matrix**

Tire	Rear Wheel		
	5 Spoke Mag	Lite Spoke	Mag Wheel
Pneumatic	20,22,24	20,22,24,26	20,22,24
Pneumatic w/ Airless Insert	20,22,24	20,22,24,26	20,22,24
Full Poly	20,22,24	20,22,24	20,22,24
High-Pressure Clincher Blk	N/A	22,24,26	N/A
Schwalbe Marathon	N/A	24,25,26	N/A

20" = 451mm 22" = 501mm 24" = 540mm 25" = 559mm 26" = 590mm

**One Arm Drive**

Available only with 24" MBL Lite spoke wheels. Std Aluminum black handrim.  
1AD option is only available on 14" - 19" frame widths.

<input type="checkbox"/> 196AL16	One Arm Drive Adj. Axle Plate	\$1,362.00
<input type="checkbox"/> 196HM1	<input type="checkbox"/> Right Hand Mount	NC
<input type="checkbox"/> 196HM2	<input type="checkbox"/> Left Hand Mount	NC

**Axle**

<input type="checkbox"/> 196AX1	Quick Release Steel	NC
<input type="checkbox"/> 196AX11	Quad Quick Release Steel <i>N/A on Spoke</i>	\$179.00
<input type="checkbox"/> 196AX45	Threaded Axle	\$148.00

**Adjustable Axle Options**

<input type="checkbox"/> 196AS4	Adjustable sleeve	NC
<input type="checkbox"/> 196AS5	Fixed Sleeve 0°	NC
<input type="checkbox"/> 196AS6	Fixed Sleeve 3°	NC

**Center of Gravity**

Measure from front of back post to center of rear axle

<input type="checkbox"/> 196CC2	1"	NC
<input type="checkbox"/> 196CC3	1.5"	NC
<input type="checkbox"/> 196CC4	2"	NC
<input type="checkbox"/> 196CC5	2.5"	NC
<input type="checkbox"/> 196CC6	3"	NC

**5. WHEEL LOCKS**

**Wheel Lock**

<input type="checkbox"/> 196WL1	High Mount Push	NC
<input type="checkbox"/> 196WL2	High Mount Pull	NC
<input type="checkbox"/> 196WL30	Compact	\$160.00
<input type="checkbox"/> 196WL24	Unilateral Wheel Lock Right	\$315.00
<input type="checkbox"/> 196WL23	Unilateral Wheel Lock Left	\$315.00
<input type="checkbox"/> 196WL10	Omit Wheel Lock	NC
<input type="checkbox"/> 196WL12	Do Not Mount	NC

**Extension Handle**

Ext. handles only available with high mount push or pull wheel locks

<input type="checkbox"/> 196EH7	Extension Handle (Pair)	\$72.00
---------------------------------	-------------------------	---------

**6. HANGERS AND FOOTRESTS**

**Front Hangers**

Footplate mounts with Extension Tube using clamping hardware

<input type="checkbox"/> 196H145	60° Swing In/Out Hanger	NC
<input type="checkbox"/> 196H139	70° Swing In/Out Hanger	NC
<input type="checkbox"/> 196H141	80° Swing In/Out Hanger	NC
<input type="checkbox"/> 196H159	Articulating Hanger <i>N/A w/ Alum Angle Adj Front Mount footplates.</i>	\$518.00
<input type="checkbox"/> 196H160	Elevating Hanger - Ext Mount <i>N/A w/ Alum Angle Adj Front Mount footplates.</i>	\$334.00
<input type="checkbox"/> 196H161	Elevating Hanger - Front Mount <i>Only available w/ front mount alum angle adj footplates and 14"-22" width</i>	\$334.00
<input type="checkbox"/> 196H144	Omit Hanger & Footplate	NC

**Articulating / Elevating Hanger Cover**

<input type="checkbox"/> 196DI47	Neoprene Cover	\$320.00
<input type="checkbox"/> 196DI48	Gel Pad Cover	\$374.00

**Extension Tubes**

See Lower Leg Length Matrix for range

<input type="checkbox"/> 196E1	Short	NC
<input type="checkbox"/> 196E2	Medium	NC
<input type="checkbox"/> 196E5	Long	NC

**Footplate**

Heel loop not included, see Footplate Option section below

<input type="checkbox"/> 196F1	Composite Footplate	NC
<input type="checkbox"/> 196F271	Composite Angle Adj-(16"-22"w)	\$173.00
<input type="checkbox"/> 196F3	Aluminum Angle Adj. Footplate	\$173.00
<input type="checkbox"/> 196F5	Aluminum <i>Only available w/Elevating Legrest</i>	NC

**Hanger & Footplate Options**

<input type="checkbox"/> 196HL12	Adjustable Heel Loop	NC
<input type="checkbox"/> 1968A02	Gel Cover For Hangers (Pair)	\$160.00

**7. BACKREST / SEATING**

**Back Type**

<input type="checkbox"/> 196BT10	Standard	NC
<input type="checkbox"/> 196BT100	Angle Adjustable Back	\$352.00

**Backrest Height**

<input type="checkbox"/> 196B5	Med (13"-17")	NC
<input type="checkbox"/> 196B9	Tall (16"-20")	NC

**Push Handles**

Integral Push Handles, Folding Height Adjustable & Omit is 8° for Standard Backrest, Straight for Angle Adj. Backrest

<input type="checkbox"/> 196PH6	Integral Push Handles	NC
<input type="checkbox"/> 196PH5	Fold-Down Push Handles <i>Only available with Adjustable Upholstery. Used in combination w/Angle Adjust only avail in 16"-20" back height</i>	\$166.00
<input type="checkbox"/> 196PH17	Folding Push Handle with Ht. Adj. Back	\$305.00
<input type="checkbox"/> 196PH2	Omit Push Handles	NC

✓	Option #	Definition	Price	✓	Option #	Definition	Price
<b>7. BACKREST / SEATING (CONT)</b>							

<b>Back Upholstery</b>			
<input type="checkbox"/>	196BU1	Standard Upholstery	NC
<input type="checkbox"/>	196BU88	Standard Upholstery EXO	\$74.00
<input type="checkbox"/>	196BU70	Backrest Tension Adj Ballistic Nylon <i>*Reference 186BU70 for billing</i>	\$187.00
<input type="checkbox"/>	196BU71	Backrest Tension Adj 3DX Vented <i>*Reference 186BU71 for billing</i>	\$464.00
<input type="checkbox"/>	196BU89	Backrest Tension Adj EXO <i>*Reference XXXBU89 for billing</i>	\$442.00
<input type="checkbox"/>	196BU14	Omit Back Upholstery	NC

<b>JAY Backs</b>			
<input type="checkbox"/>	196JB01	JAY Back (see Order forms)	Price List

<b>Military Service Patches</b>			
<input type="checkbox"/>	196VU1	US Marines Patch	NC
<input type="checkbox"/>	196VU2	US Army Patch	NC
<input type="checkbox"/>	196VU3	US Navy Patch	NC
<input type="checkbox"/>	196VU4	US Air Force Patch	NC
<input type="checkbox"/>	196VU12	US Coast Guard	NC
<input type="checkbox"/>	196VU5	Purple Heart	NC
<input type="checkbox"/>	196VU8	Purple Heart with Marines Patch	NC
<input type="checkbox"/>	196VU9	Purple Heart with Army Patch	NC
<input type="checkbox"/>	196VU10	Purple Heart with Navy Patch	NC
<input type="checkbox"/>	196VU11	Purple Heart with Air Force Patch	NC
<input type="checkbox"/>	196VU13	Purple Heart with Coast Guard	NC

<b>8. ARMREST / SIDEGUARDS</b>							
--------------------------------	--	--	--	--	--	--	--

<b>Armrest - Padded Swing Away</b>			
<i>Some limitations when used with the Angle Adjustable Back in the rear most setting</i>			
<i>Please circle a setting, will default to 11" if no selection is made</i>			
<input type="checkbox"/>	196AA18	Padded Swing-Away STD (9" to 13") Setting <input type="checkbox"/> 9" <input type="checkbox"/> 10" <input type="checkbox"/> 11" <input type="checkbox"/> 12" <input type="checkbox"/> 13"	\$274.00

<b>Armrest - Height Adjustable</b>			
<i>Flip Back HA Classic N/A w / 14"-15" depth</i>			
<input type="checkbox"/>	196AA104	Single Post Ht. Adj - Std (9.75"-14.75")	\$334.00
<input type="checkbox"/>	196AA105	Single Post Ht Adj - Low (8.25"-12.25")	\$334.00
<input type="checkbox"/>	196AA65	Flip Back Ht Adj LITE (10"-15")	\$375.00
<input type="checkbox"/>	196AA103	Dual Post Flip Back Ht Adj	\$334.00
<input type="checkbox"/>	196AA159	Flip Back Ht Adj Classic (11"-15")	\$375.00
<input type="checkbox"/>	196AA14	Omit Armrest	NC

<b>Armpad - For Height Adjustable</b>			
<input type="checkbox"/>	196AM23	Desk Length (10")	NC
<input type="checkbox"/>	196AM24	Full Length - Full (14")	NC

<b>Armrest Options</b>			
<input type="checkbox"/>	1968A01	Gel Cover For Swing Away Arms (pair)	\$160.00

<b>Side Guards</b>			
<input type="checkbox"/>	196SG33	Plastic - Standard (8.5" High)	\$185.00
<input type="checkbox"/>	196SG29	Plastic - Low (6.5" High)	\$185.00

<b>9. KEY OPTIONS</b>							
-----------------------	--	--	--	--	--	--	--

<b>Anti Tips</b>			
<input type="checkbox"/>	196AT1	Anti Tips (Pair)	\$119.00
<input type="checkbox"/>	196AT86	Anti Tips (Pair) Low Profile	\$119.00

<b>Positioning Belt</b>			
<input type="checkbox"/>	196A141	1.5" Auto Buckle w/ Clamps	\$96.00
<input type="checkbox"/>	196A804	Positioning Belt Slide Clamps	\$76.00

<b>JAY Positioning</b>			
<input type="checkbox"/>	196K1	JAY Lower Extremity S/A support <input type="checkbox"/> Left <input type="checkbox"/> Right Position	\$315.00
<input type="checkbox"/>	196K2	Omit Hanger <input type="checkbox"/> Left <input type="checkbox"/> Right	NC
<input type="checkbox"/>	196K7	Pad Size <input type="checkbox"/> 6" x 8" <input type="checkbox"/> 8" x 10"	NC

<b>Pack Accessories</b>			
<input type="checkbox"/>	196A1	Backpack - Black/Charcoal Adult	\$110.00
<input type="checkbox"/>	1961A01	Seat Pouch Black	\$82.00

<b>Spoke Guards</b>			
<input type="checkbox"/>	1964A01	Quickie Black Spoke Guard	\$155.00
<input type="checkbox"/>	1964A02	Quickie Clear Spoke Guard	\$155.00

<b>Miscellaneous</b>			
<input type="checkbox"/>	1968A03	Caddy	\$88.00
<input type="checkbox"/>	196A2100	Crutch Holder	\$113.00
<input type="checkbox"/>	196A999	Quickie Water Bottle & Holder	\$80.00

<b>10. COLOR</b>							
------------------	--	--	--	--	--	--	--

<b>Standard Frame Colors</b>			
<input type="checkbox"/>	196C140	Aztec Gold	NC
<input type="checkbox"/>	196C2	Black (Not avail as base w/Zebra)	NC
<input type="checkbox"/>	196C141	Black Cherry	NC
<input type="checkbox"/>	196C10	Black Opal	NC
<input type="checkbox"/>	196C38	Blue Opal	NC
<input type="checkbox"/>	196C36	Candy Blue	NC
<input type="checkbox"/>	196C14	Candy Purple	NC
<input type="checkbox"/>	196C35	Candy Red	NC
<input type="checkbox"/>	196C142	Electric Blue	NC
<input type="checkbox"/>	196C143	Evergreen	NC
<input type="checkbox"/>	196C44	Glossy White	NC
<input type="checkbox"/>	196C126	Glow	NC
<input type="checkbox"/>	196C144	Green Apple	NC
<input type="checkbox"/>	196C147	Hot Sparkle Pink	NC
<input type="checkbox"/>	196C148	Mauve Pink	NC
<input type="checkbox"/>	196C103	Rootbeer	NC
<input type="checkbox"/>	196C114	Sparkle Silver	NC
<input type="checkbox"/>	196C145	Sunrise Orange	NC
<input type="checkbox"/>	196C113	Titanium Color Paint	NC
<input type="checkbox"/>	196C8	Yellow	NC

<b>Matte Frame Colors</b>			
<input type="checkbox"/>	196C116	Matte Black	NC
<input type="checkbox"/>	196C151	Matte Black Cherry	NC
<input type="checkbox"/>	196C150	Matte Electric Blue	NC
<input type="checkbox"/>	196C149	Matte Evergreen	NC
<input type="checkbox"/>	196C152	Matte Purple	NC

<b>Personalized Frame Color</b>			
<input type="checkbox"/>	196C112	Desert Camouflage	\$363.00
<input type="checkbox"/>	196C43	Mossy Oak Camouflage	\$363.00
<input type="checkbox"/>	196C42	Stars & Stripes	\$363.00
<input type="checkbox"/>	196C108	Zebra (Choose above base color) Std with White unless above base paint is chosen N/A w/ Blk, Matte Black, Blk Opal, Titanium)	\$363.00

<b>Additional Options Through Built for Me</b>			
--	--	--	--

To inquire about any unique personalized Quickie products, contact: Phone: (800) 263-3390, Email: [cscanada@sunmed.com](mailto:cscanada@sunmed.com)



**APPENDIX**

**Front Seat Height Matrix Cont.**

Front Seat Height Matrix - 3/4" Stem								
Caster	3.5" FORK				5 1/4" FORK			
	STD		HEMI		STD		HEMI	
	Min	Max	Min	Max	Min	Max	Min	Max
4 Low Poly	N/A	N/A	N/A	N/A	17.0	18.0	14.0	15.0
4 x 1.25 Semi Pneu	N/A	N/A	N/A	N/A	17.0	18.0	14.5	15.5
5 Low Poly	N/A	N/A	N/A	N/A	17.5	18.5	14.5	15.5
5 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	17.5	18.5	14.5	15.5
6 Poly	N/A	N/A	N/A	N/A	17.5	18.5	15.0	16.0
6 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	17.5	18.5	15.0	16.0
8 Poly	N/A	N/A	N/A	N/A	N/A	N/A	17.0	17.0
8 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	19.5	19.5	16.5	17.0
Caster	6" FORK				7" FORK			
	STD		HEMI		STD		HEMI	
	Min	Max	Min	Max	Min	Max	Min	Max
4 Low Poly	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5
4 x 1.25 Semi Pneu	17.5	18.5	15.5	16.5	18.5	19.5	16.5	17.5
5 Low Poly	18.0	19.0	16.0	17.0	19.0	19.5	17.0	18.0
5 x 1.5 Semi Pneu	18.0	19.0	16.0	17.0	19.0	19.5	17.0	18.0
6 Poly	18.0	19.5	16.0	17.5	19.0	19.5	17.0	18.5
6 x 1.5 Semi Pneu	18.0	19.5	16.5	17.5	19.0	19.5	17.0	18.5
8 Poly	N/A	N/A	18.0	18.5	N/A	N/A	18.0	19.5
8 x 1.5 Semi Pneu	19.5	19.5	17.5	18.5	N/A	N/A	18.0	19.5

*Seat squeeze and seat depth can affect FSH availability*


Front Seat Height Matrix - 1 1/2" Stem								
Caster	3.5" FORK				5 1/4" FORK			
	STD		HEMI		STD		HEMI	
	Min	Max	Min	Max	Min	Max	Min	Max
4 Low Poly	N/A	N/A	N/A	N/A	17.5	18.5	17.5	18.5
4 x 1.25 Semi Pneu	N/A	N/A	N/A	N/A	17.5	18.5	17.5	18.5
5 Low Poly	N/A	N/A	N/A	N/A	18.0	19.0	18.0	19.0
5 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	18.0	19.0	18.0	19.0
6 Poly	N/A	N/A	N/A	N/A	18.5	19.5	18.5	19.5
6 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	18.5	19.5	18.5	19.5
8 Poly	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
8 x 1.5 Semi Pneu	N/A	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Caster	6" FORK				7" FORK			
	STD		HEMI		STD		HEMI	
	Min	Max	Min	Max	Min	Max	Min	Max
4 Low Poly	18.5	19.5	16.5	17.5	19.5	19.5	17.5	18.5
4 x 1.25 Semi Pneu	18.5	19.5	16.5	17.5	19.5	19.5	17.5	18.5
5 Low Poly	19.0	19.5	17.0	18.0	N/A	N/A	18.0	19.0
5 x 1.5 Semi Pneu	19.0	19.5	17.0	18.0	N/A	N/A	18.0	19.0
6 Poly	18.5	19.5	16.5	18.0	19.5	19.5	17.5	19.0
6 x 1.5 Semi Pneu	19.5	19.5	16.5	18.0	19.5	19.5	17.5	19.0
8 Poly	N/A	N/A	19.0	19.5	N/A	N/A	19.0	19.5
8 x 1.5 Semi Pneu	N/A	N/A	18.0	19.0	N/A	N/A	18.5	19.5

*Seat squeeze and seat depth can affect FSH availability*

**Lower Leg Length Range Matrix**

Extension Mount Hanger Leg Lengths					
Hanger	Ext Tube	Seat-to-Footer Range			
		Composite	Comp Angle Adj	Alum Angle Adj	Alum
60° Swing In/Out	Short	10.25" - 14.75"	10.25" - 14.75"	9.25" - 13.25"	N/A
	Medium	11.00" - 16.00"	11.00" - 16.00"	9.25" - 14.75"	N/A
	Long	12.50" - 16.00"	12.50" - 16.00"	10.75" - 16.25"	N/A
70° Swing In/Out	Short	11.00" - 15.75"	11.00" - 15.75"	10.00" - 14.50"	N/A
	Medium	12.00" - 16.50"	12.00" - 16.50"	10.75" - 16.50"	N/A
	Long	13.75" - 16.50"	13.75" - 16.50"	12.25" - 16.50"	N/A
80° Swing In/Out	Short	10.50" - 15.25"	10.50" - 15.25"	10.50" - 15.25"	N/A
	Medium	11.75" - 16.50"	11.75" - 16.50"	11.75" - 16.50"	N/A
	Long	13.50" - 16.50"	13.50" - 16.50"	11.75" - 16.50"	N/A

*All values are ±0.25" and based on flat seat heights and do not account for the minimum 2" ground clearance or footplate thickness. Full range not available with all seat heights (see "Min Allowable Front Seat Heights" matrix for minimum FSH information).*



**Min Allowable Front Seat Heights (for Extension Mount Hangers)**

Hanger	Ext Tube	Min Allowable FSH			
		Composite	Comp Angle Adj	Alum Angle Adj	Alum
60° Swing In/Out	Short	14.0	13.5	12.5	N/A
	Medium	14.5	14.5	12.5	N/A
	Long	16.0	16.0	13.5	N/A
70° Swing In/Out	Short	14.0	14.0	13.0	N/A
	Medium	15.0	15.0	13.5	N/A
	Long	16.5	16.5	15.0	N/A
80° Swing In/Out	Short	14.0	14.0	13.5	N/A
	Medium	14.5	14.5	14.5	N/A
	Long	16.5	16.5	15.5	N/A

**ALR & ELR Hanger Leg Lengths**

Hanger	Extension	Seat Sling to Footrest		Min Allowable FSH
		Min	Max	
ALR	Short	9 1/2	15 1/2	11 1/2
	Med	11	17 1/2	13
	Long	13	19 1/2	15
ELR Ext. Mount	Short	9	15	11
	Med	9	17	11
	Long	10	19	12
ELR Front Mount	No Ext (top hole)	4 1/2	8 1/2	6 1/2
	Short	8	11	10
	Med	10	13	12
	Long	12	15	14

Handrim	Rear Wheel		
	5 Spoke Mag	Lite Spoke	Mag Wheel
Aluminum Anodized	20,22,24	20,22,24,25,26	20,22,24
Plastic Coat	20,22,24	22,24	20,22,24
Natural Fit - Standard Grip	24	22,24,25,26	22,24
Natural Fit - No Thumb Grip	24	22,24,25,26	22,24
Natural Fit - Super Grip	24	22,24,25,26	22,24
Surge - Oval w/Gription Strip	N/A	22,24,25,26	N/A
Surge LT - Oval w/Gription Strip	N/A	22,24,25,26	N/A
Q Grip Neoprene Coated	N/A	24,25	N/A

