



AADL
ORDER FORM

February 2023

Squiggles Stander

The adjustable 3-in-1 stander designed for growing children.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		Prov	
PC		Phone	
Email / Fax			


Squiggles Stander

Weight Capacity: 48.4 lbs

✓ OPTION	DESCRIPTION	RETAIL
Standard Model		\$
Selection Required.		
Includes hip and chest support with sternum prompt; chest and hip lateral supports; size 2 knee cups, footplate and tray attachment.		
<input type="checkbox"/>	141BS34 Standard Model with Support Shell	\$2,124
<input type="checkbox"/>	141BS23 Pivot Base - Pneumatic Angle Adjustment for 3-in-1 Positioning	\$1,716
Positioning Pads		\$
Selection Required. Includes chest and hip pad covers and knee support covers.		
<input type="checkbox"/>	141ZI15 Squiggles Stander Positioning Pads - Blue	\$551
<input type="checkbox"/>	141ZI16 Squiggles Stander Positioning Pads - Blue Vinyl	\$551
<input type="checkbox"/>	141ZI17 Squiggles Stander Positioning Pads - Green	\$551
<input type="checkbox"/>	118-763 Squiggles Stander Positioning Pads - Pink	\$551
Accessories		\$
<input type="checkbox"/>	141XX31 QTY Size 1 Knee Cup	\$87
<input type="checkbox"/>	141XX32 QTY Size 3 Knee Cup	\$87
<input type="checkbox"/>	141A208 Tray - Gray	\$396
<input type="checkbox"/>	141A2165 Padded Tray Insert - Black	\$326
<input type="checkbox"/>	1418A106 Grab Rail / Toy Bar for Tray	\$321
<input type="checkbox"/>	1418A107 Small Sandals w/ Straps (Pair)	\$269
<input type="checkbox"/>	1418A108 Medium Sandals w/ Straps (Pair)	\$287
<input type="checkbox"/>	1419A24 QTY Footplate Raiser 1" (Each)	\$96
<input type="checkbox"/>	118-9500 Squiggles Stander Care Kit (Includes spare nuts, bolts, screws and washers)	\$740



Sandals



Grab Rail

Squiggles Stander

✓ OPTION	DESCRIPTION	RETAIL
Headrest		\$

Cushion selection required with hardware.

Lateral Bracket Assembly Requires Flat Headrest Support Hardware and Cushion.

<input type="checkbox"/>	141AW24	Contoured Headrest Support Hardware	\$390
<input type="checkbox"/>	141Y040	Contoured Headrest Cushion - Blue	\$191
<input type="checkbox"/>	141Y041	Contoured Headrest Cushion - Blue Vinyl	\$191
<input type="checkbox"/>	141Y042	Contoured Headrest Cushion - Green	\$191
<input type="checkbox"/>	141Y043	Contoured Headrest Cushion - Pink	\$191
<input type="checkbox"/>	141Y027	Contoured Headrest Cushion - Black	\$185
<input type="checkbox"/>	141AW22	Flat Headrest Support Hardware	\$243
<input type="checkbox"/>	141Y044	Flat Headrest Cushion - Blue	\$122
<input type="checkbox"/>	141Y045	Flat Headrest Cushion - Blue Vinyl	\$122
<input type="checkbox"/>	141Y046	Flat Headrest Cushion - Green	\$122
<input type="checkbox"/>	141Y047	Flat Headrest Cushion - Pink	\$122
<input type="checkbox"/>	141Y028	Flat Headrest Cushion - Black	\$122
<input type="checkbox"/>	141AW23	Lateral Bracket Assembly Hardware (Pair) Requires 117-697 Flat Headrest Support Hardware, Flat Headrest Cushion, and Flat Headrest Lateral Cushion	\$191
<input type="checkbox"/>	141Y048	Flat Headrest Lateral Cushions - Blue (Pair)	\$191
<input type="checkbox"/>	141Y049	Flat Headrest Lateral Cushions - Blue Vinyl (Pair)	\$191
<input type="checkbox"/>	141Y050	Flat Headrest Lateral Cushions - Green (Pair)	\$191
<input type="checkbox"/>	141Y051	Flat Headrest Lateral Cushions - Pink (Pair)	\$191
<input type="checkbox"/>	141Y052	Flat Headrest Lateral Cushions - Black (Pair)	\$191



Headrest Support



Flat Headrest



Lateral Cushions



Contoured Headrest

Technical Specifications

Frame Sizing Specification	Imperial	Metric
Age (approximate)	1 - 5 years	1 - 5 years
Maximum User Weight	48.4 lbs	22 kgs
User Standing Height	28.5" - 43.3"	730 mm - 1100 mm
Chest Width Between Laterals	6.3" - 9"	160 mm - 230 mm
Chest Support Height	21.7" - 32"	550 mm - 820 mm
Chest Support Angle Adjustment	+/- 30°	+/- 30°
Chest Support Depth Adjustment	2"	50 mm
Hip Width Between Hip Guides	6.3" - 9"	160 mm - 230 mm
Hip Support Height from Footplate to Mid Buttocks	11.8" - 20.8"	300 mm - 530 mm
Footplate Angle	10° Plantarflexion / 10° Dorsiflexion	10° Plantarflexion / 10° Dorsiflexion
Small Sandal Size	6.3" L x 3" W	175 mm L x 75 mm W
Medium Sandal Size	8.3" L x 3.2" W	210 mm L x 80 mm W
Tray Height: Prone	21.9" - 32.5"	555 mm - 825 mm
Tray Height: Supine	18.9" - 32.5"	480 mm - 825 mm
Tray Angle Adjustment	Prone 30° / Supine 40°	Prone 30° / Supine 40°
Stander Support Frame Weight	15.4 lbs	7 kgs
Stander Support Frame Dimensions	26.8" L x 15.4" W x 15.4" H	680 x 390 x 390 mm
Pivot Base Weight	12.1 lbs	5.5 kgs
Pivot Base Dimensions (Unfolded)	32.3" L x 21.7" W x 20.3" H	820 x 550 x 515 mm
Pivot Base Dimensions (Folded)	27.6" L x 21.7" W x 11.8" H	700 x 550 x 300 mm
Knee Cup Opening Width	Size 1: 2.5" Size 3: 3.5"	Size 1: 63 mm Size 3: 89 mm



Sunrise Medical Canada Inc.
237 Romina Drive, Unit 3 • Concord, ON L4K 4V3 • Canada

Customer Service: (800) 263-3390 · E-mail: cscanada@sunmed.com · www.sunrisemedical.ca

©2023 Sunrise Medical Canada Inc. MK-129857-CAN Rev (6)

