



AADL
ORDER FORM

February 2023

Squiggles Stander

The adjustable 3-in-1 stander designed for growing children.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		Prov	
PC		Phone	
Email / Fax			

Squiggles Stander

Weight Capacity: 48.4 lbs

✓ OPTION	DESCRIPTION	RETAIL
Standard Model		\$

Selection Required.

Includes hip and chest support with sternum prompt; chest and hip lateral supports; size 2 knee cups, footplate and tray attachment.

<input type="checkbox"/>	118-600	Standard Model with Support Shell	\$1,995
<input type="checkbox"/>	118-764	Pivot Base - Pneumatic Angle Adjustment for 3-in-1 Positioning	\$1,484

Positioning Pads		\$
------------------	--	----

Selection Required. Includes chest and hip pad covers and knee support covers.


<input type="checkbox"/>	118-762	Squiggles Stander Positioning Pads - Blue	\$562
<input type="checkbox"/>	118-762-03VIN	Squiggles Stander Positioning Pads - Blue Vinyl	\$562
<input type="checkbox"/>	118-760	Squiggles Stander Positioning Pads - Green	\$562
<input type="checkbox"/>	118-763	Squiggles Stander Positioning Pads - Pink	\$562

Accessories		\$
-------------	--	----

<input type="checkbox"/>	102-1640	QTY	Size 1 Knee Cup	\$111
<input type="checkbox"/>	102-3640	QTY	Size 3 Knee Cup	\$111
<input type="checkbox"/>	118-619		Tray - Gray	\$374
<input type="checkbox"/>	188-1854-05		Padded Tray Insert - Black	\$305
<input type="checkbox"/>	117-769		Grab Rail / Toy Bar for Tray	\$281
<input type="checkbox"/>	152-1600		Small Sandals w/ Straps (Pair)	\$210
<input type="checkbox"/>	152-2600		Medium Sandals w/ Straps (Pair)	\$210
<input type="checkbox"/>	120-816	QTY	Footplate Raiser 1" (Each)	\$84



Sandals



Grab Rail


Squiggles Stander

✓ OPTION	DESCRIPTION	RETAIL
Headrest		\$

Cushion selection required with hardware.

Lateral Bracket Assembly Requires Flat Headrest Support Hardware and Cushion.

<input type="checkbox"/>	120-859	Contoured Headrest Support Hardware	\$340
<input type="checkbox"/>	137-800-03	Contoured Headrest Cushion - Blue	\$167
<input type="checkbox"/>	137-800-03VIN	Contoured Headrest Cushion - Blue Vinyl	\$167
<input type="checkbox"/>	137-800-01	Contoured Headrest Cushion - Green	\$167
<input type="checkbox"/>	137-800-04	Contoured Headrest Cushion - Pink	\$167
<input type="checkbox"/>	137-800-05	Contoured Headrest Cushion - Black	\$167
<input type="checkbox"/>	117-697	Flat Headrest Support Hardware	\$228
<input type="checkbox"/>	117-850-03	Flat Headrest Cushion - Blue	\$114
<input type="checkbox"/>	117-850-03VIN	Flat Headrest Cushion - Blue Vinyl	\$114
<input type="checkbox"/>	117-850-01	Flat Headrest Cushion - Green	\$114
<input type="checkbox"/>	117-850-04	Flat Headrest Cushion - Pink	\$114
<input type="checkbox"/>	117-850-05	Flat Headrest Cushion - Black	\$114
<input type="checkbox"/>	118-667	Lateral Bracket Assembly Hardware (Pair) Requires 117-697 Flat Headrest Support Hardware, Flat Headrest Cushion, and Flat Headrest Lateral Cushion	\$162
<input type="checkbox"/>	137-807-03	Flat Headrest Lateral Cushions - Blue (Pair)	\$298
<input type="checkbox"/>	137-807-03VIN	Flat Headrest Lateral Cushions - Blue Vinyl (Pair)	\$298
<input type="checkbox"/>	137-807-01	Flat Headrest Lateral Cushions - Green (Pair)	\$298
<input type="checkbox"/>	137-807-04	Flat Headrest Lateral Cushions - Pink (Pair)	\$298
<input type="checkbox"/>	137-807-05	Flat Headrest Lateral Cushions - Black (Pair)	\$298



Headrest Support



Flat Headrest



Lateral Cushions



Contoured Headrest

Technical Specifications			
Frame Sizing Specification	Imperial		Metric
Age (approximate)	1 - 5 years		1 - 5 years
Maximum User Weight	48.4 lbs		22 kgs
User Standing Height	28.5" - 43.3"		730 mm - 1100 mm
Chest Width Between Laterals	6.3" - 9"		160 mm - 230 mm
Chest Support Height	21.7" - 32"		550 mm - 820 mm
Chest Support Angle Adjustment	+/- 30°		+/- 30°
Chest Support Depth Adjustment	2"		50 mm
Hip Width Between Hip Guides	6.3" - 9"		160 mm - 230 mm
Hip Support Height from Footplate to Mid Buttocks	11.8" - 20.8"		300 mm - 530 mm
Footplate Angle	10° Plantarflexion / 10° Dorsiflexion		10° Plantarflexion / 10° Dorsiflexion
Small Sandal Size	6.3" L x 3" W		175 mm L x 75 mm W
Medium Sandal Size	8.3" L x 3.2" W		210 mm L x 80 mm W
Tray Height: Prone	21.9" - 32.5"		555 mm - 825 mm
Tray Height: Supine	18.9" - 32.5"		480 mm - 825 mm
Tray Angle Adjustment	Prone 30° / Supine 40°		Prone 30° / Supine 40°
Stander Support Frame Weight	15.4 lbs		7 kgs
Stander Support Frame Dimensions	26.8" L x 15.4" W x 15.4" H		680 x 390 x 390 mm
Pivot Base Weight	12.1 lbs		5.5 kgs
Pivot Base Dimensions (Unfolded)	32.3" L x 21.7" W x 20.3" H		820 x 550 x 515 mm
Pivot Base Dimensions (Folded)	27.6" L x 21.7" W x 11.8" H		700 x 550 x 300 mm
Knee Cup Opening Width	Size 1: 2.5"	Size 3: 3.5"	Size 1: 63 mm Size 3: 89 mm



Sunrise Medical Canada Inc.

237 Romina Drive, Unit 3 • Concord, ON L4K 4V3 • Canada

Customer Service: (800) 263-3390 · E-mail: cscanada@sunmed.com · www.sunrisemedical.ca

©2023 Sunrise Medical Canada Inc. MK-129857-AADL Rev (4)

