

ORDER FORM




Q200R with Rehab Seat

EIPW35

Customer		Ship To	
Account #		Name	
Date		Address	
PO#		City	
Buyer		Province	Postal Code
Mark For		Phone	

Client Measurements

a) Top of head to seat	in
b) Top of shoulder to seat	in
c) Underarm to seat	in
d) Elbow to seat	in
e) Top of leg to seat	in
f) Shoulder width	in
g) Chest width	in
i) Hip width	in
j) Actual sitting depth	in
m) Foot length	in
n) Knee to foot	in
o) Seat width	in



Notes

Chair Specifications

Weights starting at: (with Batteries)	213 lbs
Power Seating	22° Tilt
Seat Type	Rehab Seat w/ Fold Down Back and 12° Angle Adjustment
Seat Size Range	Width 15" to 21" Depth 16" to 20"
Seat-to-Floor Height Range	17" to 19"

Mid-Wheel Drive Base Specifications

Drive Wheel Position	Rear-wheel drive
User Weight Capacity	300 lbs.
Base Size	Width 23", Length 42"
Minimum Turning Radius	18.5" (excludes footplates)
Drive Control	R-net PM 90A / VR-2 90A
Driving Range	Up to 15.5 miles
Drive Motors	2-pole, 6 mph / 10 kph
Drive Wheels	12" w/ Flat Free Inserts 14" w/ Pneumatic Inserts 14" w/ Flat Free Inserts
Front Casters	9" w/ Pneumatic Inserts 9" w/ Flat Free Inserts
Charger	8A Dual Mode (Compact), 10A
Battery Size	Group 22
Color	7 color accents
Ground Clearance	3"

BASE			
OPTION		DESCRIPTION	RETAIL
Standard Duty Rear Wheel Drive Base <i>End Item Number EIPW35 (Rehab Seat)</i>			\$
350PC302		Quickie Q200 R NPC No Power Option Included in Base Pricing : Q200 R Base, VR-2 90amp Electronics, 4-wheel Suspension System, 12" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 10kph / 6mph, Weight Capacity = 300 lbs., STFH = 17".	\$7,467.00
350PC303		Quickie Q200 R SPC Single Power Option Included in Base Pricing : Q200 R Base, VR-2 90amp Electronics, 4-wheel Suspension System, 12" Drive Wheels with Dual Density Inserts, 9" Flat Free Casters, ISO Occupied Tie-Down Brackets and 8amp Charger. Max Speed = 10kph / 6mph, Weight Capacity = 300 lbs., STFH = 17".	\$7,467.00
User Weight			\$
350WG23		Up to 300 lbs.	NC
Motor Package			\$
350PF107		2-pole 6 mph / 10 kph	NC
Drive Wheels			\$
350DW77		Drive Wheel, 12", Solid, Knobby, Black	NC
350DW78		Drive Wheel, 14", Pneumatic, Knobby, Black	NC
350DW79		Drive Wheel, 14", Solid, Knobby, Black	NC
Caster Wheels			\$
350CW121		Caster Wheels 9", Pneumatic, Black	NC
350CW122		Caster Wheels 9", Solid, Black	NC
Color Accents			\$
350EZ20		Apple Green	NC
350EZ22		Black	NC
350EZ21		Blue	NC
350EZ23		Orange	NC
350EZ24		Red	NC
350EZ19		White	NC
350EZ25		Silver	NC
ISO Transit Tie Downs			\$
350T050		ISO Transit Package	NC
LED Lights			\$
350LI04		LED Lights and Indicators	\$1,207.00
350LI02		Omit LED Lights	NC
350LI08		LED Fender Lights	\$636.00

BASE			
OPTION		DESCRIPTION	RETAIL
Seat Heights <i>(Seat height is measured from front of seat pan to floor with 0° frame angle.)</i>			\$
350HT55		Seat Height 17" (N/A with Tilt)	NC
350HT56		Seat Height 18"	NC
350HT57		Seat Height 19" (N/A with Tilt)	NC
Batteries and Chargers			\$
Group 22 Battery			
350AG31		GP22 Flat Top Pair w/Straps <i>(N/A w/ Tilt)</i>	\$668.00
350AG35		GP22/50Ah Flat Top Pair w/Straps	\$668.00
Chargers			
350BJ30		10A Dual Mode (Compact)	\$160.00
350BJ19		8A Dual Mode (Compact)	NC
SEATING			
OPTION		DESCRIPTION	RETAIL
Power Seating Options			\$
Seat Tilt			
350PA56		Posterior 22° Tilt	\$3,582.00
Pre Set Seat Frame Angle <i>(No Power Seating Options, Increases Seat to Floor Height)</i>			\$
350FA10		0° <i>(Adds 0" to STFH)</i>	NC
350FA3		3° <i>(Adds 1" to STFH)</i>	NC
350FA4		6° <i>(Adds 2" to STFH) (N/A w/ 19" STFH)</i>	NC
Backrest Frame Type			\$
350BF27		Tool Angle Adjustable Back, 90° to 102°, Stroller Handle	NC
350BF30		Tool Angle Adjustable Back, 90° to 102°, Trailing Push Handles	NC
Backrest Frame Option			\$
350BF31		Fold Down, Quick Release Backrest	NC
Back Cane Height			\$
350B266		Tall, 23"	NC
Back Cushion Type			\$
350BT148		Lite Sling Back, with Tension Adjustable Straps	NC
350BT135		JAY® J3® (Section A)	NC
350BT133		Omit Back Cushion	NC

SEATING			
OPTION		DESCRIPTION	RETAIL
A - JAY® J3 Back			\$
JAY® J3 Back Contour			
350JP01		J3 PA Back for 14" to 20" Wide	\$899.00
350JP11		J3 PL Back for 14" to 20" Wide	\$971.00
350JP13		J3 PD Back for 14" to 20" Wide <i>(18" and 21" shell height only)</i>	\$1,052.00
JAY® J3 Shell Height			
350JR14		16.5" Tall Shell <i>(Upper Thoracic Short)</i>	NC
350JR15		18" Tall Shell <i>(Upper Thoracic Med)</i>	NC
350JR16		19.5" Tall Shell <i>(Upper Thoracic Tall)</i>	NC
350JR17		21" Tall Shell <i>(Shoulder Height)</i>	NC
350JR27		22.5" Tall Shell <i>(Shoulder Height)</i>	NC
350JR28		24" Tall Shell <i>(Shoulder Height)</i>	NC
JAY® Anterior Thoracic Support			
350K10		Medium <i>(ATS8302)</i>	\$317.00
350K11		Large <i>(ATS8303)</i>	\$317.00
JAY® Lateral Thoracic Swing-Away Supports			
350JW35		Lateral Thoracic (LxH) 100 mm x 127 mm / 4" x 5" Pad - Right	\$71.00
350JW38		Lateral Thoracic (LxH) 100 mm x 127 mm / 4" x 5" Pad - Left	\$71.00
350JW36		Lateral Thoracic (LxH) 100 mm x 152 mm / 4" x 6" Pad - Right	\$71.00
350JW39		Lateral Thoracic (LxH) 100 mm x 152 mm / 4" x 6" Pad - Left	\$71.00
350JW37		Lateral Thoracic (LxH) 152 mm x 100 mm / 6" x 4" Pad - Right	\$71.00
350JW40		Lateral Thoracic (LxH) 152 mm x 100 mm / 6" x 4" Pad - Left	\$71.00
Lateral Thoracic Swing-Away Hardware			
350JW41		Lateral Thoracic Swing Away Bracket - Right	\$317.00
350JW42		Lateral Thoracic Swing Away Bracket - Left	\$317.00
Seat Type			\$
350ST80		JAY® <i>(Comes with Solid Seat Pan) (See Section B)</i>	NC
350ST34		SEDEO Cushion <i>(Sizes 16", 18", & 20" Wide) (Comes with Solid Seat Pan)</i>	\$421.00
350ST85		Solid Seat Pan	NC
350ST84		Seat Sling	NC

SEATING			
OPTION		DESCRIPTION	RETAIL
Headrests (Requires Pad, Upholstery and Mount)			\$
Headrest Pad			
350HE36		SEDEO Pad Small With Cover (L x H) 240 x 150 mm / 10" x 6" <i>(Only Available with SEDEO Multi-Position Headrest Mount and SEDEO Upholstery)</i>	\$224.00
350HE38		Whitmyer® Plush Headrest Pad 203 mm / 8"	\$371.00
350HE39		Whitmyer® Plush Headrest Pad 254 mm / 10"	\$371.00
350HE44		Whitmyer® Plush Headrest Pad 356 mm / 14"	\$371.00
350HE46		Whitmyer® Adjust-A-Plush Headrest Pad Medium <i>(Adjustable from 7" to 16" Wide x 4" Tall)</i>	\$438.00
350HE48		OMIT Headrest Pad	NC
Headrest Mount			
350MH32		SEDEO Multi-Position Headrest Mount with Quick Release <i>(Only available with SEDEO Pad)</i>	\$549.00
350MH33		Whitmyer® Linx Headrest Mount with Quick Release	\$243.00
350MH34		Whitmyer® Linx ² Headrest Mount with Quick Release	\$294.00
350MH35		Whitmyer® AXYS Headrest Mount with Quick Release	\$243.00
350MH17		OMIT Headrest Mount	NC
Seat Size			\$
Seat Pan Width (Seat width is determined by inside of arm pads.)			
350R78		15" <i>(16" seat pan width)</i>	NC
350R53		16" <i>(16" seat pan width)</i>	NC
350R54		17" <i>(16" seat pan width)</i>	NC
350R55		18" <i>(19" seat pan width)</i>	NC
350R56		19" <i>(19" seat pan width)</i>	NC
350R57		20" <i>(19" seat pan width)</i>	NC
350R58		21" <i>(19" seat pan width)</i>	NC
Seat Pan Depth (Front of Seat Pan to Back Structure.)			
350RF51		16"	NC
350RF45		17"	NC
350RF46		18"	NC
350RF47		19"	NC
350RF48		20"	NC
B - JAY® Seat Cushion			\$
350CU59		JAY® Ion® <i>(Sizes 16" - 21" Wide)</i>	\$498.00
350CU61		JAY® Union® <i>(Sizes 16" - 21" Wide)</i>	\$753.00
350CU64		JAY® Fusion® <i>(Sizes 16" - 21" Wide)</i>	1,030.00
350CU104		JAY® Fusion® with Cryo® Technology	\$1,135.00

POSITIONING SUPPORTS & ACCESSORIES

OPTION	DESCRIPTION	RETAIL
Armrest Structure		\$
350EG42	Single Post Removable <i>(Only available with Single Power, SPC)</i>	NC
350EG43	Single Post Removable	NC
Side Guards (Pair)		\$
350SG49	Side Guards <i>(Depth Adjustable)</i>	\$351.00
Arm Pads		\$
350AM15	2" x 14" Armpad	NC
Swing-Away Legrest (Pair)		\$
350H177	SEDEO PRO Legrests, Swing Away, 70° <i>(14" to 21" Leg Length)</i>	NC
350H178	SEDEO PRO Legrests, Manual Elevating, Swing Away <i>(Includes Padded Calf Pads) (14".5 to 21.5" Leg Length)</i>	\$374.00
350H179	SEDEO PRO Legrests, Power Elevating, Swing Away <i>(Includes Padded Calf Pads) (15.5" to 22.5" Leg Length)</i>	\$3,384.00
350H180	SEDEO PRO Hanger Receivers Only	NC
Swing Away Legrest Calf Supports		
350DI42	Calf Support, Depth Adjustable, Padded <i>(6.5" to 8" Wide, Depending On Chair Width, x 9.5" Tall)</i>	\$504.00
Footplates		
350F370	Angle Adjustable Flip Up Footplates	\$280.00
Accessories		
350DI44	Heel Loops	NC
350KG4	Knee Pads, Integral Foam <i>(2.5" Wide x 4.5" Tall)</i>	\$278.00
SEDEO Lower Extremity Swing Away Support		
350VD01	Lower Extremity Support, Padded, Swing Away, SEDEO, Right	\$504.00
350VD02	Lower Extremity Support, Padded, Swing Away, SEDEO, Left	\$504.00
Pelvic Positioning Belts		
350A837	Positioning Belt, Auto Style Buckle, 50 mm 1600 mm / 2"	NC
350A839	Positioning Belt, Aircraft Buckle, 50 mm x 1300 mm / 2"	NC
350A840	Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2"	\$101.00
Accessories		\$
Backpack		
350A1	Black Utility	\$110.00
Medial Thigh Support (N/A with Sling Seat)		
3508A80	Medial Thigh Support	\$267.00
3508A81	Pommel Swing-Away Hardware	\$384.00
Misc. Accessories		
3508A52	Ram Mount Phone	\$76.00
3508A53	Ram Mount Tablet	\$303.00
3508A54	Ram Mount Cup Holder	\$76.00
350A611	Crutch Holder	\$133.00
3508A58	Rear View Mirror	\$76.00
3509A10	Utility Hooks	\$57.00

JOYSTICK / ELECTRONICS

OPTION	DESCRIPTION	RETAIL
VR-2 Joystick \$		
350J127	VR-2 Joystick	NC
R-net Joystick (Remote) \$		
350J121	Standard LED Joystick (w/o Screen Display)	NC
350J138	Color Joystick (2.3" Color LCD Screen Display, CJSM1)	NC
350J118	High Resolution Color Joystick, (CJSM2+BT) 3.5" Color LCD Display, Paddle Switches, Programmable Shortcut Keys, ICONs, Stereo Jack, IR and dual (IOS/Android) Bluetooth control	\$1,245.00
350J119	Omit Remote Control	NC
Joystick Mounting \$		
350JM34	Fixed on Armrest	NC
350JM35	Parallel Swing-Away	\$440.00
Joystick Mount (Location) \$		
350HM16	Right Hand Armrest	NC
350HM17	Left Hand Armrest	NC
Display Language Selection R-net \$		
350L1	English	NC
350L3	French	NC
350L7	Spanish	NC
Joystick Handle (Optional) N/A with HD Joystick \$		
350JH13	Ball Handle	\$108.00
350JH15	T-Handle	\$108.00
Electronic Programmings \$		
350PR10	R-net PC Programmer w/ Dongle	\$696.00
R-net Power Modules \$		
350ER12	Power Module, VR-2 90A, Non-Expandable	NC
350ER13	Power Module, R-net, 90A, Expandable	NC
350EX15	Expandable 90A Controller - R-net	NC
350EX16	Harness for 90A Expandable Controller	NC

JOYSTICK / ELECTRONICS

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

Seating Control

\$

Powered Seating Control

350E040	Omit Thru-Drive <i>(Requires Switch Control and Standard LED Joystick Requires Ctrl+5) (N/A w/ VR-2 electronics)</i>	NC
350E037	Thru-drive Control for 1 Actuator	\$1,422.00
350E038	Thru-drive Control for 2+ Actuators	\$2,860.00

Switch Control *(Available with Omit Thru-Drive only) (Requires use of Omni2 or Joystick Jacks)*

350E033	Dual Toggle for 1 Actuator	NC
350E034	2.5" Single Button Switch for 1 Actuator <i>(Button switches are not mounted)</i>	NC
350E035	Switch-It Mini 1 button for 1 Actuator	NC
350E036	Switch-It Mini 2 button for 1 Actuator	NC
350E041	Ctrl+5 Switch Box for 2 Actuators or more	NC

Additional Power Seating Control

350QG11	Ctrl+5 Switch Box with Free Assignable Functions	\$1,062.00
---------	--	------------

Dual Toggle Mounting Location

350TN1	Toggle Mount, Right	NC
350TN2	Toggle Mount, Left	NC

Ctrl+5 Mounting Location

350EY19	Right Armrest	NC
350EY20	Left Armrest	NC
350EY21	Back Cane, Right	NC
350EY24	Back Cane, Left	NC

Ctrl+5 Switch Type

350EY22	Push Buttons	NC
350EY23	Toggle Switches	\$351.00

Attendant Controls

\$

350A05	Dual Control Module for Attendant	\$1,207.00
350A06	R-net Color LCD (CJSM1), Attendant Control Only	\$1,207.00

Attendant Mounting Location

\$

350DN5	Right Back Cane	NC
350DN6	Left Back Cane	NC

R-net Programming

\$

Assignable Buttons

Assignable Button feature allows for customization of joystick buttons and/or Ctrl+5 Switches allowing access to various wheelchair and seating functions. When enabled, chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

350PM12	Include Assignable Buttons	NC
---------	----------------------------	----

JOYSTICK / ELECTRONICS

OPTION	DESCRIPTION	RETAIL
--------	-------------	--------

Specialty Input Controls

\$

R-net Omni2

350ON14	R-net Omni2 Only Input Control Module and 3.5" LCD Display w/ Built-In Infrared & Bluetooth Choose Mount Below	\$2,349.00
350ON18	R-net Omni2 Bundled Input Control Module and 3.5" LCD Display w/ Built-In Infrared & Bluetooth Bundled Omni2 Available w/ Switch-It Input Devices	\$769.00

Omni2 Mount (If used w/ Standard Joystick, it must be mounted opposite side)

350DS14	Link-it with Parallel Swing Away Omni2 Mount w/ Rear Mounted Interface, Right	\$688.00
350DS15	Link-it with Parallel Swing Away Omni2 Mount w/ Rear Mounted Interface, Left	\$688.00

Proportional Hand Input Devices

350EU14	Switch-It® MicroGuide	\$4,574.00
350EU15	Switch-It® MicroPilot	\$4,574.00

Hand Input Device Mounting Position (Required)

350EM05	Link-It to the Fixed Joystick Mount	\$410.00
350EM06	Link-It to the Parallel Swing- Away Joystick Mount	\$410.00
350EM08	Omit Mount	NC

Hand Input Device Mounting Platforms (Optional)

350MG16	SunDisk Hand Tray <i>(only available with MicroSeries joysticks)</i>	\$410.00
350MG15	Bullet Hand Tray	\$410.00

Hand Input Device Mounting Location (Required)

350IM09	Right	NC
350IM10	Left	NC

Proportional Head Array Input Devices

350EU25	Switch-It® DualPro Proportional Head Array - Facial Pads	\$6,953.00
350EU26	Switch-It® DualPro Proportional Head Array - Spot Pads	\$6,953.00

Head Array Mounting (Required)

350IJ13	Whitmyer® Linx Headrest Mount with Quick Release	\$243.00
350IJ14	Whitmyer® Linx ² Headrest Mount with Quick Release	\$294.00
350IJ15	Whitmyer® AXYS Headrest Mount with Quick Release	\$243.00



Q200 R with Rehab Seat

\$

Base Price

Total of additional equipment

Packaging and transport

Total Price

Date

Signature

Month

Day

Year



Customer Service: (800) 263-3390
Fax: (800) 561-5834
E-mail: cscanada@sunmed.com
www.sunrisemedical.ca
Specifications are subject to change without notice.

© 2024 Sunrise Medical Canada Inc.
237 Romina Drive, Unit 3
Concord, ON
L4K 4V3

MK-129793 CAN Rev. 2